



HEALTH & BEAUTY

ENGLISH

2023



www.promed.de



Certification DIN EN ISO 13485

Valid from 1. Januar 2023

TABLE OF CONTENTS

Electric files

Disinfection	5
Files for professional nail design	7
UV lamps & dust extraction	17
Files for home use and travel	20
Bits	90

Health control

Masks	4
Air purifier	6
Blood pressure monitors	24
Thermometers	27
Pulse oximeter	29
Blood glucose monitor	31

Therapy

Inhalation	27
TENS/EMS	32

Wellness

Massage	34
Heat	37
Ambient air / Aroma diffuser.....	43
Fragrance oils	46

Shoes, soles and accessories

Promed comfortable shoes introduction	48
Overview of shoes and indication	52
Promed Size charts.....	54
Promed "Pedibelle" comfortable shoes	55
Promed post-op shoes	57
Huber soles by Promed	66
Darco	68
Darco Size charts.....	88

Accessories

Bits amm.	90
Arm Cuffs + Inhalation	106
Electrodes + TENS/EMS	108
Medisana Accessories	110

Advertising materials	111
------------------------------------	-----

Glossar	112
----------------------	-----

Imprint	115
----------------------	-----



Art. 333001 | EAN 40 43641 33023 1

MNS-100

Particle-Filtering Half Mask
Standard EN 149:2001 +A1:2009;
Regulation (EU) 2016/425

promed

Also approved for studios, salons,
hairstylists, pedicurists...

Efficient three-layer filtration system:
spunbond, meltbond, non-woven



Protection class
FFP2 NR



3D design for com-
fortable breathing



3-layer
filtration system



One size
fits all

- Protection against substances potentially harmful to human health, such as fine particulate matter, liquid aerosols and biological substances (bacteria, viruses, fungi). Effective prevention of allergens and germs. Protects you and your environment.
- High-grade filtration of 95% of all particles in the air
- Skin-friendly material with high wearing comfort
- Tested and CE-certified according to standard EN 149:2001 +A1:2009; Regulation (EU) 2016/425 on personal protective equipment

- Optimum fit: V-shape with nose contour
- Hygienic individual packaging prevents the mask from being contaminated before it is used. Single foiled in PE-bags.
- One size fits all. Masks for onetime use. Do not wash.
- Content: 10 pieces, incl. hook and instructions
- Also available in larger quantities on request – for information, see: info@promed.de

TECHNICAL DATA

Type of mask: FFP2 NR, category III - Product size (folded): 106.7 x 155 mm - Product weight: approx. 8 g - Shelf life: 3 years, CE 2834

MNS-100 Black

Particle-Filtering Half Mask
Standard EN 149:2001 +A1:2009;
Regulation (EU) 2016/425

Efficient three-layer filtration system:
spunbond, meltbond, non-woven

Also approved for studios, salons,
hairdressers, pedicurists...



Protection class
FFP2 NR



3D design for com-
fortable breathing



3-layer
filtration system



One size
fits all

- Protection against substances potentially harmful to human health, such as fine particulate matter, liquid aerosols and biological substances (bacteria, viruses, fungi). Effective prevention of allergens and germs. Protects you and your environment.
- High-grade filtration of 95 % of all particles in the air
- Skin-friendly material with high wearing comfort
- Tested and CE-certified according to standard EN 149:2001 +A1:2009; Regulation (EU) 2016/425 on personal protective equipment

- Optimum fit: Vshape with nose contour
- Hygienic individual packaging prevents the mask from being contaminated before it is used. Single foiled in PE-bags.
- One size fits all. Masks for onetime use. Do not wash.
- Content: 10 pieces, incl. hook and instructions
- Also available in larger quantities on request – for information, see: info@promed.de

TECHNICAL DATA

Type of mask: FFP2 NR, category III - Product size (folded): 106.7 x 155 mm - Product weight: approx. 8 g - Shelf life: 3 years, CE 2834



Disinfection for instruments (concentrate)

Des sum exerci musapid eritatum,
quatem. Maionecatur atqui nos
dolorerit, sincita tureptia t

Coming soon:
New product in 2023

- Highly effective, aldehyde-free concentrate for the disinfection/cleaning of instruments and sensitive rotary instruments
- Atiosanditas vit offic to blam voloraest faccaec taquam facit, omnine lat omniat atur sitibus doluptur, omni ad quati o
- doloruntem excepratem faceptaquist aut lacius most arum landelestiam quosaperit laborem que que custin cus a dolorum aut abo. Bus pos aligendita aut quid qui nest qui

- a. Aces eat volest, esti sed moloreroriae nonsequiate est, iusciatis abo. Ro excereria cus sequatur aut accum lam
- Laboruntia serum eos mo blandip issimus etur?
- Sedis que quo eos verum facese velitincta nonsece rferrume
- Vitasped ut et optaspedit ratiati amende nonserc hitibus el molestrum deligna tianihit quo doluptas voluptatibus sitaec

Art. 334000 | EAN 40 43641 33026 2

AC-4000

Air purifier with Wlan function

Cleaning power CADR
of up to 420 m³/h

Effective against viruses & bacteria,
all kinds of dust, mould, pollen, mites,
smoke/odours



Possible to operate
by smartphone app



Remote control included




especially quiet

- High filter capacity due to combination filter (HEPA filter H13, activated carbon, photocatalytic, cold catalytic)
- 3-level aluminium pre-filter (washable)
- UV disinfection by means of ultraviolet light
- 3-stage colour display AQI for air quality (Air Quality Indicator), integrated particle sensor
- Cleaning power CADR (Clean Air Delivery Rate) of up to 420 m³/h
- Timer function of 1 to 9 hours
- Temperature indicator and humidity indicator
- 3 power levels of the fan as well as sleep mode, auto mode, manual mode
- Suitable for rooms up to 70 m²
- Automatic reminder to replace the filter
- Easy to use thanks to large touch display, child safety lock

Promed AC-4000 – the ideal solution for clean air in closed rooms.

This high-quality air purifier with combination filter effectively filters pollutants such as viruses, bacteria, all kinds of dust (including fine dust), mould, pollen, mites, smoke and odours. Up to 99.95% of the frequently occurring microorganisms are removed from the air.



- ✓ Filters particles especially effectively with an efficiency of up to 99.95%
- ✓ Improved air purification by means of an ion generator: generates 5 million anions/cm³;
- ✓ Improves air quality and promotes health, particles are safely removed

TECHNICAL DATA

Application area:	approx. up to 70 m ²	Power supply:	220 - 240 VAC, 50 Hz Rated output: 85 W
Cleaning power CADR:	380 - 420 m ³ /h (Clean Air Delivery Rate)	Noise level during operation:	35 dB (A) - 55 dB (A)
Anions Concentration:	generates 5 million/cm ³ (improved air quality)	Dimensions in cm:	approx. 41.5 x 25.5 x 69.5
UV lamp:	220 V, 4 W, 254 nm, UV-C, socket T5/65	Weight:	approx. 8 kg

SCOPE OF DELIVERY

- 1x air purifier AC-4000
- 1x aluminium pre-filter and combination filter
- 1x UVC disinfecting tubular lamp
- 1x remote control
- 2x 1.5 V batteries AAA
- 1x instructions for use

[Art. 334001]



[Art. 334002]

Replacement accessories
can be ordered:
Filter set and UVC lamp





Advantages of Promed files

We offer the latest technology and best service from one source.

Premium quality for long-lasting products. Our highly qualified staff is always there for you in case of questions.



Professional equipment for skin, nail and foot care

As a leading manufacturer of technically matured and proven professional cutters, Promed develop, research and design its devices in close cooperation with users and medical technology engineers.

Our more than 30 years of experience allows us to add meaningful developments for you, which facilitate you to work faster and comfortable on the customers.

Promed grew as the first manufacturer for technical innovations such as ASC (Automatic Speed Control) or the proven OPC (Overload Protection Control).

Promed devices are regularly tested by recognized testing institutes and undergo a highest German and global safety standards. Our strengths couch to outstanding technical innovations unique in after-sales service (only in Germany).

We are represented at the following beauty trade fairs in Germany:





Technology that makes the difference

APE – Automatic Power Electronics

For powerful, effective operation. The intelligent electronics supply the motor with more electricity even at low rpm.

ASC – Automatic Speed Control

The digital automatic speed control registers the loss of speed when pressure is exerted during filing and automatically compensates for this, allowing the unit to maintain its set power. You always have the full power at your disposal, letting you work even more quickly.

OPC – Overload Protection Control

The unique overload protection control automatically registers an overload of the respective motor and thus ensures that the device is not damaged. It is available for handpieces and, if applicable, for suction turbines.

Overload protection control for the handpiece motor:

In the event of a blockage or an overload of the handpiece, the handpiece and, if applicable, the suction turbine are switched off.

Overload protection control for the suction turbine:

If the suction turbine is not cooled adequately and the air temperature in the suction turbine housing exceeds a critical value, the suction turbine and the handpiece are switched off.

Safety Stop

If the device develops excessive speed and pressure during grinding, the electronics will respond by switching off the device. In this case, you must set the rotation switch of the device to zero (middle position) and start the device again.

Surfaces

Nail polish remover resistant foils with a stylish silver metallic look.

Performance

Better performance is achieved for greater engine power by intelligent, microprocessor-controlled electronics.

All-in-One control

Large, handy and easy to use, this allows you to select speed and rotation direction simultaneously.

Handunits

A new development made from a lightweight-construction aluminum-stainless steel combination with precision bearing. Extra-powerful micro motors. Superb ergonomics and vibration-free operation. Non-metal hand units are made from special synthetic material.

Promed – benefit yourself from our Know-How.



manicure



ASC (Automatic Speed Control)



Safety Stop



pedicure



OPC (Overload Protection Control)



Cold-light lamp/LED Lamp



artificial nails



APE (Automatic Power Electronics)



Storage battery operation

the file - 5040-SX2

Professional heavy duty spraying unit
for foot care in wet technique

NEW: with a color display!

All types of filing work in the foot care.
Nail and cornea treatment.



- 4 memory keys for the programming of individually set speeds
- Powerful marked micro pumps for a finest spray dust (adjustable)
- Ergonomic handpiece with an integrated coolant spray device
- Automatic nozzle cleaning during switching on and off
- Automatic pump cleaning

- Overload protection for the motor of the handpiece engine
- Stand-by button for all functions
- Stand-by button on the handpiece
- Water tank (410ml) – with a slightly refillable plug-in system
- Optional: foot pedal (on/off) and electronic foot controller (see accessories p. 104)

TECHNICAL DATA

Power supply: 230/115 VAC, 50/60 Hz, 105 W · 207 x 163 x 245 mm (L x B x H) · Pump capacity: 550 mbar; 8 ml/min · Weight: Handpiece 112 g, total unit approx. 3.2 kg



NEW FEATURES OF THE DEVICE

(as compared to the predecessor)

① Automatic nozzle cleaning:

When turning the device on, the device only blows air through the nozzles first to remove any dirt. When turning it off, the spray pump continues running for several seconds to clean the nozzle and/or the hose system.



② Spray setting

The file is equipped with a high-performance micropump generating finest spray. It offers the possibility of individually adjusting the spray pattern of the spray to the spray liquid used.



Scale for setting the spray pattern

③ Abrasive product recommendation:

The abrasive product recommendation of the device indicates which maximum speed is to be worked with for which milling head diameter.



④ Automatic pump cleaning:



If work has been carried out for a period of more than one hour with the spray function switched off, the device automatically activates the pump cleaning. This process is shown on the display and usually takes 2 - 3 seconds.



the file - 4030-SX2

Love your health!

Modern file with vacuum device

Particularly low dust exposure thanks to direct suction on the abrasives. Perfect for nail design and foot care.



product film:



- Membrane keyboard (soft keys) with digital display
- Continuously adjustable up to 30,000 rpm (clockwise/anti-clockwise rotation)
- 3 memory buttons for individuelle speed
- Display indicating the required replacement of the filter bag
- Temperature-related safety switch-off
- Powerful suction capacity (continuously adjustable 350 watts turbine), very low-noise
- Flexible suction hose inkl. suspension

- Featherweight ergonomic handpiece (95 g) with standby switch
- Incl. 5+1 filter bags (14x18 cm), easy to change, tear-resistant
- Incl. handpiece holder
- Optional: foot pedal (see accessories on p. 104)

TECHNICAL DATA

Dimensions: (W x H x D) 225 x 170 x 230 mm - Weight: approx. 2.6 kg - Supply voltage: 230 VAC/50 Hz, 350 VA (also available as 115 VAC/60 Hz, 350 VA) - Turbine: 350 watts - Noise emission: 62 dB A, speed: up to 30,000 rpm

Vacuum performance and number of revolutions



Keyboard view and filter replacement display



WORKING IN A HEALTH-CONSCIOUS MANNER

In order not to risk your health, the fine dust is exhausted directly at the handpiece during work. The occurrence of dust in the studio is reduced to a minimum. The practical guiding of the suction hose is ensured by our unique suspension at the device.



ACCESSORIES



Filter bag 5x (white) + 1x in the device



Handpiece holder



Cleaning tool for the collet



Tube suspension

OPTIONAL ACCESSORIES



Foot switch (on/off) or foot pedal (swell)

the file - 2020

Latest file generation for continuous operation

For professional treatment in
continuous operation.



- Rotational speed is infinitely variable up to 20.000 rpm
- Ergonomic, powerful, vibration-less stainless steel handpiece (160g) especially for pedicure
- Professional rotating collet
- 5 LEDs for rotational speed indicator
- 1 LED for foot pedal indicator
- Safety Stop; ASC (Automatic speed control)
- All-in-One regulator



- Rotary knob with bright, sparkling gemstone (diamond shape)
- Foot pedal optionally available (on/off) (see accessories on p. 104)

TECHNICAL DATA

Rotational speed: infinitely variable up to 20,000 rpm · Power supply: 100 - 240 VAC, 50/60 Hz, 800 mA;
Sec. 24 VDC, 1.5 A · Weight: 450 g · Weight of handpiece: 160 g · Dimensions: 130 x 200 x 100 mm (D x W x H)

the file - 1030

Latest generation of file
for professional nail design

For a professional treatment of natural
and artificial nails as well as nail design
in continuous operation



product film:



- Rotational speed is infinitely variable up to 30,000 rpm
- Ergonomic, powerful, vibration-free stainless steel handpiece (160g)
- Professional rotating collet
- 5 LEDs for rotational speed indicator
- 1 LED for foot pedal indicator
- Safety Stop; ASC (Automatic speed control)



- All-in-One regulator
- Rotary knob with bright, sparkling gemstone (diamond shape)
- Foot pedal optionally available (on/off) (see p. 104)

TECHNICAL DATA

Rotational speed: infinitely variable up to 30,000 rpm · Power supply: 100 - 240 VAC, 50/60 Hz, 800 mA;
Sec. 24 VDC, 1.5 A · Weight: 450 g · Weight of handpiece: 160 g · Dimensions: 130 x 200 x 100 mm (D x W x H)

the file - 625

For using manicures, simple
naildesign and even gentle pedicures



- Rotational speed is infinitely variable up to 25,000 rpm
- Light ergonomic, powerful, vibration-free stainless steel hand piece (160 g)
- Professional rotating collet
- 5 LEDs for rotational speed indicator
- All-in-One regulator (clockwise and counterclockwise, On/Off)

- Safety Stop
- Rotary knob with bright, sparkling gemstone (diamond shape)


TECHNICAL DATA

Power supply: 100 - 240 VAC, 50 / 60 Hz, 1200 mA; Sec. 24 VDC, 1.5 A - Energy consumption with zero load: 0.3 W - Rotation: 25 000 rpm - Measurements: approx. 220 x 120 x 100 mm (L x B x H) - Weight: Handpiece 95 g

the file - 620

The entry-level device
for ambitious nail design

Is also well suited for finger and
foot nail, no matter natural, gel or acrylic.
You may also remove hard skin and
callus easily.



- Space-saving, super lightweight body
- Handpiece in purple-glitter look (75 g)
- Automatic quick release clamp for standard bits (Shaft ø 2.332 - 2.350 mm)
- Support for the handpiece on the unit
- All-in-One regulator (20,000 rpm.; clockwise/anti-clockwise; On/Off)

- Rotary knob with bright, sparkling gemstone (diamond shape)
- Safety-Stop (Stops at overloading); Speed display
- Lightest device on the market!

TECHNICAL DATA

Power supply: 100 - 240 VAC, 50 / 60 Hz, 800 mA; Sec. 24 VDC, 1.5 A - Energy consumption with zero load: 0.27 W - Rotation: 20,000 rpm. - Measurements: approx. 110 x 90 x 80 mm (L x B x H) - Weight: Handpiece 75 g



Art. 330210 | EAN 40 43641 33004 0

UC-50 Ultra sonic device

promed

Use for: genuine jewellery,
glasses, denture, instruments,
grinding tools, CD/DVDs etc.

Possible use with tap water



product film:



Figure is only used as an example

- Cleans instruments, grinding tools, genuine jewellery, watches, CD/DVDs, third set of teeth etc. with ultrasonic high-frequency vibrations
- The UC-50 switches off automatically after the timer expires
- Very stable and quiet

- Illuminated LCD display; timer function (5 timer settings)
- Suitable for tap water
- Suitable for Promed Pure-ID (improves cleaning and disinfecting)
- Accessories: basket, bit holder, CD holder, watch holder

nailfan

Unit for removing harmful fine dust

Tested dust protection class "L" -
DIN EN ISO 60335-2-69



product film:



No cold fingers due to the lateral removal of cold air

- The Promed nailfan stands out due to its tested dust protection class "L", quality, easy handling and special performance
- Ergonomic design for non-tiring work – for you and your clients
- Removable antistatic padded arm rest
- No cold fingers despite powerful suction
- Very large dust collection area of up to 20 cm
- Low-noise and low power consumption

- Large (11 x 19 cm), easy to change fine dust particle filter (washable!)
- No filter bags required
- Separate dust drawer; Antistatic housing
- Incl. second reusable filter, towel for catching dust
- Accessories: see p. 105


TECHNICAL DATA

Main unit dimensions: (B x H x T) 308 x 140 x 340 mm - Weight: ca. 1.8 kg - Power supply: 230 VAC, 50 Hz - Fan: Power rating: 18/22 W, Noise level: 46 dbA - Max. air flow 180 m³/h

nailfan mini

Dust extraction

Ideal for the studio and at home



- Very lightweight, easy-to-use device
- Extraction adjustable between 2 levels (low-noise)
- High extraction capacity on the entire support surface
- Filter (reusable) easy to change and clean
- Ergonomic design with customer friendly support surface which customers can rest their hands on

- Small, handy and portable – also ideal for on the move
- Including extra filter
- Accessories: see p. 105

TECHNICAL DATA

Voltage supply: 100 – 240 VAC, 50/60 Hz 0.6 A, sec: 12 V / 1.5 A - Fan: Level I: 3,000 rpm / 24 W; Level II: 4,500 rpm / 36 W - Dimensions: (W x H x D) 240 x 185 x 60 mm - Weight of the device: approx. 570 g

Art. 223005 | EAN 40 43641 22004 4

promed

UVL-54 „All-In“

UV LED Light Curing Device



UV-LED lamp for hardening
all materials



- Suitable for fingernails and toenails to cure gel and acrylic
- Dual power system with 36 LEDs for even curing
- 4 different timer settings (10s/30s/60s/99s)
- Anti-heat function with 99s timer
- High definition LEDs: Service life approx 50,000 hours
- Removable base (magnetic)



- Automatic sensor system (on/off – integrated light barrier) till 120s
- Digital display

TECHNICAL DATA

Voltage supply: 100 VAC - 240 VAC; 50/60 Hz, 0.3 A, Sec. 12 VDC, 3 A - Wavelength of UV LED (Å): 365 nm - 405 nm - LED service life: approx. 50 000 h - Dimensions (LxWxH) in mm: 221 x 205 x 92

Art. 330040 | EAN 40 43641 33019 4 (white)

Art. 330045 | EAN 40 43641 33020 0 (pink)

promed

UVL-36 S

UV light for curing gel and acrylic

Incl. 4x UV lights
(easily interchangeable)




- Professional nail design for the home
- Ideal for curing gel and acrylic
- 360 degree reflection for even curing
- Easy operation thanks to a 4-stage time (90s/120s/30min/60min)
- Removable base for easy replacement of the light

- 4 slip-proof rubber feet for stable standing
- See page 105 for accessories

TECHNICAL DATA

Dimensions: (L x W x H) approx. 253 x 221 x 113 mm - Weight: approx. 1.04 kg - Power supply: 230 VAC, 50 Hz, 36 W - Operating periods: 90 s / 120 s / 30 min / 60 min - Light socket G23 - Wavelength of the lamps: λ = 365 nm



LTM-30 LED table magnifying lamp – for working comfortably



- Compact LED table magnifying lamp for your workplace, your workshop, for reading, and much more
- 30 individual LED bulbs ensure excellent light output
- 175% magnification while working (3 dpt)
- High quality 4 inch lens with protective cover
- Easy adjustment of the lens position and height
- Swivelling head and arm joints (multidimensional)
- Ergonomic design and stable base

- Extremely long-life LEDs (high-quality long-life)
- LEDs are much more robust than conventional neon tubes
- Cool shadowless lighting for bright, comfortable working
- Lens size: 4-Inch lens (Ø 10.2 cm)
- Colour temperature: 6,200 K (neutral white)
- Dimensions (ØxVxH): 16 x 18 x 39 mm

LTL-749 LED table lamp – for maximum flexibility



- Classy and stylish design
- Very light and even illumination
- No tired eyes, relaxed working
- Ideal for your workplace, on any table, in the bedroom or at school
- 3-level dimming of brightness at the touch of a button
- Multidimensional 3-axis adjustment
- LED bar adjustable by 360°

- High-quality long-life LED tube: appr. 40,000 hours
- Energy-saving compared to conventional lamps
- LED lamp with integrated glare protection for your eyes
- Frame made of high-quality aluminium
- Anti-slip base
- Dimensions of the LED bar: approx. 33 cm long
- Brightness: approx. 5,000 K (colour temperature)

Art. 200630 | EAN 40 43641 20033 6

620 deluxe

(manicure/pedicure)



Incl. 43 accessories

The lightest device on the store!

„The profi class
quality at home“



- From nail cosmetics and design to the treatment of thickened skin
- Chuck (with new engine) with purple-glitter look (75g light)
- Separate power supply for maximum power
- All-in-One controller; Overload protection
- Rotary knob with bright, sparkling gemstone (diamond shape)

Art. 114010 | EAN 40 43641 11006 2

Pedimed (horny skin remover)

Professional pedicure set for silky soft feet




- Efficiently and gently remove excess skin easily.
- Easily and simply moves excess horny skin on your feet and heel from the first application
- Soft and efficient treatment of your skin and nails
- Safe and easy operation due to its ergonomic design
- Slim, lightweight file fits comfortably in your hand
- Powerful 20,000 revolutions per minute (switchable clockwise/anti-clockwise rotation)
- High-quality, durable abrasive wheels with a standard shaft dental diameter
- Also ideal when on travel and on the go
- Incl. 30 abrasive belts, 5 abrasive wheels in a designer box, designer case

Ultra pro S (manicure/pedicure)

Incl. 42 accessories

"Elegant top class"



- Satisfies with its comfortable operation
- The grinding force can be continuously adjusted up to 13,000 rpm
- The professional right/left rotation makes it easier to treat difficult areas

- An all-rounder, from satisfying nail care and design to the removal of calluses
- Including: storage box

Pedisenso Duo (manicure/pedicure)

Incl. 44 accessories

„The file for diabetics“

product film:



- Injury-free treatment thanks to the „diabetic safety chuck“
- Automatic stop when too much pressure is applied
- Established, powerful technology using epicyclical gearing (for calluses removal)

- Grinding force is continually adjustable up to 4,000 rpm
- Accessories: see p. 102 for collets
- Including: storage box

Art. 173010 | EAN 40 43641 17001 1

promed


Sensitive

(manicure/pedicure)

Incl. 42 accessories

„Top class quality
with LED-light“

product film:



- Cold light illumination (with LED-lighting), Security Stop
- Comfortable, dust proof soft touch operation buttons control clockwise/counterclockwise rotation
- Continuously adjustable up to 12,000 rpm
- Power safety transformer, intensity display with 5 LEDs
- Incl. storage box and microfiber duster

Art. 100030 | EAN 40 43641 10005 6

promed

Emotion

(manicure/pedicure)

Incl. 41 accessories

„All-rounder with LED-light“

inkl. storage case



with cold light illumination



- For the simple and safe care of nails and problem zones, such as thickened skin on the ball of the foot and the heel
- Continuously adjustable up to 12,000 rpm
- Safety stop prevents possible injury if too much pressure is applied
- Comfortable operation using buttons for grinding force and with right/left rotation

Feeling (manicure/pedicure)

Incl. 6 accessories

„Accu powered“



inkl. bit holder



- Trims and shapes the nails, repairs cracked (artificial) nails, removes cuticles and thickened areas of skin very easy
- New chuck to change the bits faster, optimized rotation-balance

- 2 speed levels (7,000rpm, 11,000rpm)
- Maximal freedom of movement without cable
- Very powerful motor (new), for use with accu





Health and blood pressure



Blood pressure monitors support better preventative health care and optimise your health.

The first step to better health: monitoring

Health checks increasingly become the focus of many people – on the one hand, due to the better health education of the population and, on the other, because of the progress in medical research.

Regularly measure and check precisely your blood pressure and heart rate. Monitor your health and identify irregularities in due time.



PBM-4.2

Upper arm blood pressure monitor

incl. cuff (22-42 cm),
storage bag and 4 batteries
WHO Classification
 WHO=World Health Organization

- Precise measurement at your upper arm
- 2x90 memory spaces (2 users)
- Display screen: systolic, diastolic, heart rate, date, time
- Arrhythmia indicator
- Average value display of the last 3 measurements
- Assessment of the measured values according to the WHO



Cuffs in different sizes see p. 106

- Large digital LCD screen
- Automatic switch-off function after 3 minutes
- Includes cuff size 22 - 42 cm,
- Includes 4xAA batteries (1.5 V)
- Measuring range: – 0-299 mmHg;
– Pulse: 40-199 beats/min

PBM-3.5

Upper arm blood pressure monitor

incl. cuff, size M (22 - 32 cm),
storage bag and batteries
WHO Classification
 WHO=World Health Organization

- Precise blood pressure measurement on the upper arm
- 90 memories for each of 2 users
- Display readouts: systolic, diastolic, pulse, date, time
- Arrhythmia display
- Average value display of the last 3 measurements
- Large digital LCD display: approx. 72x53 mm



Cuffs in different sizes see p. 106

- Automatically switches off after 3 min
- Operation with 4 x AA batteries or with 6 Volt AC Adaptor
- Including: Storage bag and batteries
- Measuring range: – Blood pressure: 0-299 mmHg
– Pulse: 40-199 beats/min

HGP-40

Wrist blood pressure monitor



WHO Classification
WHO=World Health Organization

- 60 memory spaces for 1 user
- Display: systole, diastole, heart rate, date, time
- Arrhythmia indicator
- Automatic switch-off function after 3 minutes of inactivity
- Classification of the measured values according to the classification system of the World Health Organisation

- Average display
- Including: 2 x AAA (1,5 VDC) batteries,
- Storage box
- Cuff for 135-215 mm wrist circumference
- Measuring range:
 - Blood pressure: 40-230 mmHg
 - Heart rate: 40-199 beats/min

HGP-30

Wrist blood pressure monitor



WHO Classification
WHO=World Health Organization

- 120 memory spaces each for 2 users
- Display: systole, diastole, heart rate, date, time
- Arrhythmia indicator
- Calculation of the average of the last 3 saved value
- Classification of the measured values according to the classification system of the World Health Organization

- Battery level indicator
- Includes: 2xAAA (1.5V) batteries, storage box
- Cuff for 135 - 210 mm wrist size
- Measurement range:
 - Blood pressure: 30-280 mmHg
 - Heart rate: 40-200 beats/min

INH-2.1 Ultrasound inhalator

Gentle treatment in the case of respiratory diseases

Suitable for allergy sufferers: 100 % latex-free

Including handy starter set for on the move




- High efficiency due to innovative micro membrane nebulisation (mesh technology) with tiny droplets
- Short inhalation time due to very high nebulisation performance (0.32 ml/min)
- Particularly efficient due to the high respirable aerosol proportion
- Inhalation with mask or mouthpiece – also easy to use for children
- Low-noise nebulisation

- Small and handy (portable) – ideal for on the move
- Filling volume content: max. 8 ml
- Nebulisation performance: approx. 0.32 ml/min
- Particle size: 2.1 μm
- Size: approx. 42x45x93 mm
- Including: storage bag, mouthpiece, mask for children and adults, connecting piece, batteries
- Optionally available: power supply unit (Art. 422006)
- Accessories: see p. 107

IRT-100

Infrared thermometer



Precise measuring within seconds



Normal temperature display



"Fever alarm" through colour indication

- Non-contact infrared clinical thermometer
- Hygienic and safe measuring without skin contact
- Within seconds: precise measuring of the body, room and surface temperature
- Visual fever alarm; 30 memory spaces
- 3 different display colours for normal temperature, slight fever and fever

- Display in °C and F
- Automatic switch-off (after 30sec)
- Measuring range for body: 32.0-42.9° C
- Measuring range for surfaces/objects: 0-100° C
- Including 2xbatteries

Art. 432020 | EAN 40 43641 40078 1

IRT-80

Infrared thermometer

promed



Precise measuring within seconds



- Non-contact infrared fever thermometer
- Hygienic and safe measuring
- Measurement of body and object temperature within seconds
- Fever alarm (visual and acoustic)
- 35 memory spaces
- 2 different display colours

- For normal temperature and fever
- Display in °C and °F
- Automatic switch-off (after 10 sec.)
- Measuring range for body: 32.0-42.9 °C
- Measuring range for surfaces/objects: 0-100 °C
- Includes storage bag and 2 batteries



Normal temperature display



"Fever alarm" through colour indication

Art. 431005 | EAN 40 43641 40076 7

PFT-3.7

Digital clinical thermometer

promed

Including:
storage box, battery

- Precise measurement of your body temperature: orally, axillary, rectally
- Acoustic signal when ready for measurement
- Acoustic temperature alarm
- Last measured value is automatically saved

- Waterproof
- Automatic switch-off after 10 minutes
- Digital LCD display
- Easy battery change
- Measuring range: 32.0-42.9 °C

Art. 481010 | EAN 40 43641 40082 8

PM-200 Pro

Fingertip pulse oximeter

Pulse oximeter for measuring blood oxygen levels (SpO₂%) and heart rate (pulse frequency)



- Display in column or wave format
- 6 different displays to choose from
- 4 adjustable brightness levels
- Two-colour OLED display
- Battery level indicator

NEW

promed



Display in portrait and landscape

Includes carrying strap

- Automatic switch-off after 8 seconds
- Includes carrying strap and batteries
- Dimensions: approx. 59.4x33x31.2 mm
- SpO₂: Measurement range: 70%-100%;
- PR: Measurement range: 30 beats/min – 254 beats/min

Art. M79456 | EAN 40 15588 79456 8

medisana.

PM 100 connect

Pulse oximeter



- Data transfer via Bluetooth® Smart (4.0) to the VitaDock+ app for iOS and Android*
- Pulse oximeter to measure the saturation of oxygen in blood (SpO₂%) and the heart rate (pulse)
- View measurements on OLED display
- Display brightness adjustable in 10 levels
- Simple „one touch“ operation – ideal for on the go
- Oxygen saturation (SpO₂) and pulse rate can be displayed as bars or waves
- 6 different views
- Automatic switch-off after 8 seconds
- Low battery display
- Incl. batteries (2x 1,5V AAA)

Art. M79455 | EAN 40 15588 79455 1

medisana.

PM 100

Pulse oximeter



- Pulse oximeter to measure the saturation of oxygen in blood (SpO₂%) and the heart rate (pulse)
- Oxygen saturation (SpO₂) and pulse rate can be displayed as bars or waves
- View measurements on a OLED display
- Viewable from 6 different angles
- Display brightness adjustable over 10 levels
- Simple „one touch“ operation – also ideal for travelling
- Automatic switch-off after 8 seconds
- Battery level indicator
- Including: Batteries, carrying strap

Art. M79030 | EAN 40 15588 79030 0 | mg/dL
 Art. M79034 | EAN 40 15588 79034 8 | mmol/L

medisana.

Meditouch 2 USB

Blood glucose monitor

Measuring time within 5 seconds

Including „starter set“:
 10 test strips, lancing device, 10 lancet,
 AST cap, control solution, 5 alcohol pads,
 pouch, quick reference guide,
 diabetic diary, batteries



starter set

- Precise measurement accuracy according to the new ISO 15197: 2013 thanks to new test strip technology (GDH-FAD)
- Underfilling detection
- USB port for data transfer to VitaDock Online and to the VitaDock app for iOS (7.0) and Android (4.3)
- Pre and post meal marker
- Acoustic test reminders
- Calculates average values for 7/14/30/90 days
- 480 memory slots
- Automatic ketone warning
- Starter set incl. lancing device, AST cap, 10 lancets, 10 test strips, control solution, bag
- 0,6µl blood
- Extra large, easy-to-read display
- Accessories, e.g. test strips see p. 110

Art. M79048 | EAN 40 15588 79048 5

medisana.

Meditouch 2 connect DUAL USB

Blood glucose monitor

Including a „starter set“:
 10 test stripes, lancing device,
 10 lancets, AST cap, control solution,
 bag, user instruction manual, batteries

Can only be ordered in Germany!



starter set

- Precise measuring accuracy corresponds to the new ISO 15197:2013 thanks to the new test stripe technology (GDA-FAD)
- Underfilling detection
- Acoustic test reminder
- Measuring unit individually adjustable: mg/dL and mmol/L
- Memory for 480 measuring results
- Bluetooth® Smart (4.0) data transfer to VitaDock + App for iOS (from 7.0) and Android (from 4.3)
- Starter set incl. lancing device, AST cap, 10 lancets, 10 test stripes, control solution, bag
- Accessories, e.g. test stripes, see p. 110

Art. 352430 | EAN 40 43641 35034 5

IT-6 Pelvic floor trainer

Therapeutic device for training and strengthening the pelvic floor



Incl. 4 electrodes, 1x vaginal probe, Batteries and practical storage box



2
YEARS
WARRANTY

- Used to prevent or reduce bladder weakness
- Treatment of incontinence after giving birth
- Training the sphincter muscles – for treating stress/urge/mixed/faecal incontinence
- Using the electrical stimulation, you regain a feel for the pelvic floor and thus learn to effectively contract the pelvic muscles again

- The device fulfils all your need through simultaneous treatment of different regions using 4 electrode connections; thanks to the dual outputs, two probes can be used at the same time and controlled separately
- For anal and vaginal stimulation, can improve the sex life
- See p. 108 for available accessories

Art. 352420 | EAN 40 43641 35030 7

EMT-6 Combi-unit

Digital, electrical pain therapy and muscle stimulation unit

Can also be used in sport



product film:



6x
EMS
TENS
5x

2
YEARS
WARRANTY

- Unit for effective pain and muscle therapy
- It is predestined for the treatment of chronic and acute pain conditions, mobility impairments and for rehabilitation, for the preservation, building up and training of muscles and muscle groups (e. g during sport)
- It can be prescribed by the doctor (cash refund)



- Two-channel unit, biphasic-asymmetric square pulse current, DC-neutral
- Package contents: combined unit for TENS & EMS in a practical storage case, 4 electrodes (40x40mm), 2 electrode cables, batteries, beltclip
- Different electrodes see accessories p. 108

EMT-4

Combi-unit for electrical nerve and muscle stimulation, wellness



- 11 preset programs (5 Tens and 6 EMS programs)
- Intensity and power-on protection
- Combi-unit for Electrical Nerve and Muscle Stimulation, pain therapy and wellness
- Clock function in standby mode
- EMT-4 is able of order (cash refund)
- Scope of delivery: Promed EMT-4, 2x self-adhesive electrodes (40x40 mm), unit pedestal, batteries, 2x electrodes cable, manual instruction, keeping case, waist bag
- Different electrodes see accessories p. 108

Tens 1000 s

Analog pain therapy device



- Apart adjustable frequency and impulse width, three amplitude modulations standing for the choice
- Timer: continuously, 15 and 30 min adjustable.
- The Promed tens 1000 s is able of order (cash refund).
- Scope of delivery: Promed Tens 1000s, 4 self-adhesive electrodes (40x40 mm), 9-Volt battery, 2 contact cable, instructions, a practical plastic suitcase
- Different electrodes see accessories p. 108



Art. 621010 | EAN 40 43641 60006 8


promed

ITM Pro

Hand massage device



2 attachments each for finger massage
and for acupressure



incl. 2 different attachments



- Powerful tapping massage action
- Smoothly adjustable control of massage intensity
- Connectible red light and warming function
- Various massage attachments: Acupressure and pressure point massage

- Ergonomic shape makes it easy to use
- Non-slip handle
- Extra-long wire
- 230VAC / 50Hz / 35 W


MSA-900

Shiatsu massage seat cover



Washable, removable neck rest cover

Height-adjustable neck massage



Height-adjustable



- Targeted, intensive high-definition massage effect for a relaxed sense of well-being
- Ergonomic fit of the massage seat cover
- Individually adjustable neck massage height
- 3 massage zones that can be chosen: entire, upper or lower back
- Separate vibration massage in the seat with 3 selectable intensity levels

- Additionally connectible red light and warming function
- Automatic switch-off after 15 minutes
- Easy-to-clean, durable material
- Washable, removable neck rest cover
- Practical hook and loop fastener for fixation to the chair back
- Power supply: 220-240 VAC / 50-60 Hz / 48 W

FB-100 Foot spa



(1)



(2)



(3)

Including 3 pedicure top parts



exchangeable



- Relaxing bubble and vibration massage
- Pleasant red light field
- Keep-warm function
- Massage rollers for foot reflex zone massage
- Massage-supporting footbed
- Dry massage function

- 3 exchangeable pedicure top parts: (1) brush, (2) massage top part, (3) cornea remover
- Slip-proof rubber feet
- Comfortable size (suitable up to shoe size 46)
- Power supply: 220-240 VAC / 50-60 Hz, 80 W

Art. 340010 | EAN 40 43641 34001 8

promed

Intimate douche „Orchidee“

Device for feminine hygiene



- Ergonomic design with easy handling
- Stainless steel
- 7 jets at the head of the rod
- Incl. 2x pressure regulators and 1x storage box

- Connection thread: G1/2 inch
- Autoclavable and dishwasher-safe
- Made in Germany
- 10 years warranty

Art. 340020 | EAN 40 43641 34002 5

promed

Intimate douche „Rose“

Device for feminine hygiene



- Ergonomic design with easy handling
- Rosfree, stainless steel
- 10 jets at the head of the rod
- Incl. line pressure regulator and manual

- Autoclavable and dishwasher-safe
- 2 years warranty
- Connection thread: G1/2 inch
- Weight: approx. 205 g

Excellent heating products.
Excellent quality.



KHP-2.3 G

Heating blanket "Sherpa"

NEW

promed



- 10 temperature levels
- 9 different adjustable periods of time (1-9 hours)
- Quick, efficient heating up with 150 watts (Turbo Heat)
- Safety system with overheat protection
- Automatic switch-off after selected period of time set
- Comfortably cosy warmth on the entire surface
- Practical control element

- Fluffy, soft outer material, two-sided
- Can be washed at 40°C thanks to removable switch
- Material: 100% polyester flannel
- Size: approx. 180x130 cm
- Power: 150 watts

KHP-2.3

Cosy heating blanket

promed



- 10 temperature levels
- 9 different adjustable times (1-9 hours)
- Fast, powerful heating with 150 watts (turbo heat)
- Safety system with overheating protection
- Automatic switch-off after selected time setting
- Pleasant soft heat on the entire surface

- Practical control element
- Fluffy, soft outer material
- Washable at 40°C thanks to removable switch
- Material: 100% polyester flannel
- Size: approx. 180x130 cm
- Output: 150 watts

NRP-2.0

Shoulder and neck heating pad



- 6 temperature levels
- Automatic switch-off after 90 minutes
- Heats up quickly with 100watts (Turbo Heat)
- Heats into the neck area
- Safety system with overheating protection
- Fluffy-soft upper material which is comfortable to wear

- Easy to use, ergonomic operation/switch
- Washable at 30°C
- Material: 100% polyester
- Size: approx. 52x56 cm
- Power: 100 watts

NRP-5.4

Neck and Back Heating Pad



- 4 temperature levels
- Automatic switch-off after 90 minutes
- Heats up quickly with 100watts (Turbo Heat)
- Heats into the neck area
- Safety system with overheating protection

- Easy to use, ergonomic operation/switch
- Fluffy-soft upper material which is pleasant to wear
- Material: 100% polyester
- Size: approx. 60x90 cm
- Output: 100 watts

NRP-2.4

Neck and back heating pad



- 10 temperature levels
- 9 different adjustable times (10-90 min)
- Fast heating with 100 watts (turbo heat)
- Automatic switch-off
- Safety system with overheating protection
- Pleasantly soft, removable cover
- Machine washable at 40°C

- Particularly soft plush fleece cushion
- Perfect support due to back & neck-friendly fit
- Practical Velcro fastener for an individual adjustment of pad
- Material: 95% plush fleece, 5% nylon
- Size: approx. 39x63 cm
- Output: 100 watts

HGP-1.7

Heating belt



- 10 temperature levels
- 9 different adjustable times (10-90 min)
- Fast heating with 100 watts (turbo heat)
- Safety system with overheating protection
- Perfect support due to adjustable Velcro fastener
- Ergonomic operation/switch
- Suitable for an abdominal girth up to approx. 125 cm

- Variable wearing comfort: around the abdomen or diagonal across the back
- Practical Velcro fastener for an individual adjustment of pad
- Washable at 40°C thanks to removable switch
- Size: approx. 58x30 cm
- Material: 90% plush fleece, 10% nylon
- Output: 100 watts

HKP-1.2 H

Cosy heating pad



- 3 temperature levels
- Automatic switch-off after 90 minutes
- Cover and pad are washable at 40°C thanks to removable switch
- Safety system with overheating protection
- Easy handling, ergonomic operation/switch
- Fast heating with 100watts (turbo heat)

- Includes pad cover (also usable without cover)
- Fluffy soft outer material and supply cover
- Material: 100% polyester
- Size: approx. 40x30 cm
- Weight: approx. 0.4 kg
- Output: 100 watts
- Power supply: 230 VAC / 50 Hz

HKP-1.6 XL

Heating pad



- 6 temperature levels
- Automatic switch-off after 90 minutes
- Heats up quickly with 100watts (Turbo Heat)
- Safety system with overheating protection
- Easy to use, ergonomic operation/switch
- Includes a 2-sided heading pad cover (can also be used without cover)

- Fluffy soft outer material and cuddly cover
- Cover washable at 30°C
- Material: 100% polyester
- Size: approx. 46x35 cm
- Output: 100watts

KFS-100

Foot warmer

"No more cold feet"



- 6 temperature levels
- Turbo Heat with powerful 100 watt
- Automatic switch-off after 90 minutes
- Safety system with overheating protection
- Ergonomic control element with illuminated LED control display

- Removable, fluffy-soft inner lining
- Inner lining can be machine-washed at 30°C
- Particularly soft, breathable fabric core
- Large foot space (suitable up to shoe size 46)
- Material: microfibre (foot warmer), flannel (inner lining)
- Size: approx. 30x30x23 cm

Art. 452210 | EAN 40 43641 40058 3

UBS-2.4

Underblanket



- 10 temperature settings
- Short warming up time thanks to 120 watts (turbo heat)
- Safety system with overheating protection
- Automatic switch-off after selected time (1-9 hours)
- Pleasantly gentle and evenly distributed heat for a relaxing sleep
- Convenient control element with display

- Cosy breathable material
- Machine washable at 40 °C thanks to removable switch
- Material: 100 % cotton (top), 100 % Polyester (bottom)
- Size: approx. 150 x 70 cm
- Output: 120 watts



Art. 651025 | EAN 40 43641 60014 3

AL-1 300 WS

Aroma diffuser with 7-fold LED colour change

promed

Wellness



Water tank and nebulizer zone are off antibacterial material – a possible bacteria emissions are effectively reduced

Very large water tank XXL: 1.3 litre

Automatic switch-off when the tank is empty



- Aroma diffuser in noble wood look, timeless design suitable for every manner of room design; ideal for aromatic oils
- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere
- Using ultrasound technology to create a fine mist of fragrance that spreads throughout the entire room
- Highly effective due to micro-fine atomisation

- Running time approx. 10-14 hours, strong nebulisation
- Suitable for all large and medium-sized rooms
- Wellness light with 7-fold colour change
- Low-noise and energy-saving
- Including: Power adapter
- 100-240 VAC / 50-60 Hz
- 24VDC / 500mA / 12 W
- Flavours, see p. 46/47

promed

AL-240 FL

Aroma diffuser with flame light

NEW

promed

LED wellness light with flame effect

Dimensions: approx. 210 x 135 x 81 mm



- Aroma diffuser with flame light for aromatic oils
- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere
- Uses ultrasonic technology to create a fine mist of fragrance that spreads throughout the room
- High efficiency due to micro-fine atomisation
- Automatic switch-off when the tank is empty

- Low-noise and energy-saving
- Including: mains adapter and water dosing cup
- Large 240ml water tank
- 100-240 VAC / 50-60Hz
- 24VDC / 500mA / 12W
- Flavours, see p.46/47



AL-400 CL

Aroma diffuser with integrated clock and 7-fold LED colour change

NEW

promed

Decorative & elegant wood look features

Dimensions: approx. 200 x 135 x 143 mm



- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere – Ideally suited for scented oils
- Using ultrasonic technology, creates a fine mist of fragrance that spreads throughout the room
- Wellness light with automatic colour change in 7 colours (colours: can be individually set and dimmed)
- Interval function (30 sec.) of the scented mist

- Highly effective due to micro-fine atomisation
- Automatic switch-off when the tank is empty
- Low-noise and energy-saving
- Including: power adapter and water measuring cup
- Large 400ml water tank
- 100-240 VAC / 50-60 Hz, 300mA
- 24VDC / 500mA / 12W
- Flavours, see p.46/47



AL-150 WS

Aroma diffuser with 7-fold LED colour change

Aroma Diffuser in high-quality wood look

Automatic switch-off when the tank is empty



- Humidifies and freshens the air in the room for a soothing and therapeutic atmosphere
- Using ultrasound technology to create a fine mist of fragrance that spreads throughout the entire room
- Highly effective due to micro-fine atomisation
- Large and efficient 150ml water tank

- Wellness light with 7-fold colour change
- Low-noise and energy-saving
- Including: Power adapter
- 100-240 VAC / 50-60 Hz
- 24VDC / 500 mA / 7 W
- Flavours, see p. 46/47

AL-305

Aroma diffuser with 5-fold LED colour change

300ml water tank – runtime of approx. 5-9 hrs

Automatic switch-off when the tank is empty



- Aroma diffuser for fragrance oils
- Using ultrasound technology, it creates a fine mist of fragrance that spreads throughout your entire room
- High efficiency due to micro-fine atomisation
- Automatic switch-off when the tank is empty
- Timer to program the operating time: 2 and 4 hours
- Wellness light with 5-fold colour change

- Low-noise and energy-saving
- Including: refillable reservoir, power adapter
- 300ml water tank; runtime: approx. 5-9 hrs
- 100-240 VAC / 50-60 Hz
- 24VDC / 500 mA / 12 W
- Flavours, see p. 46/47

Art. see below



Fragrance oil / Aroma essences

(100% essential oil / composition / perfume oil)



Fragrance oil Relax (10 ml)

Art. 262010

EAN 40 43641 26005 7



Fragrance oil Be Happy (10 ml)

Art. 262011

EAN 40 43641 26006 4



Aroma essences Lavender (10 ml)

Art. 262012

EAN 40 43641 26007 1



Aroma essences Lemongrass (10 ml)

Art. 262013

EAN 40 43641 26008 8



Aroma essences Orange (10 ml)

Art. 262015

EAN 40 43641 26010 1



Nature-identical oil, cinnamon (10 ml)

Art. 262017

EAN 40 43641 26012 5



Aroma essences Eucalyptus (10 ml)

Art. 262018

EAN 40 43641 26013 2



Promed fragrances are made in Germany – from the best natural products processed to obtain 100 % essential oils – and spread a pleasant and sensual scent. Indulge yourself.

Fragrance oil / Aroma essences

(100 % essential oil / composition / perfume oil)



Aroma essence cinnamon-orange (10 ml)

Art. 262019

EAN 4043641 26015 6



Aroma essence lemon (10 ml)

Art. 262020

EAN 4043641 26016 3



Perfume oil peach (10 ml)

Art. 262021

EAN 4043641 26017 0



Perfume oil apple-cinnamon (10 ml)

Art. 262022

EAN 4043641 26018 7



Perfume oil wild cherry (10 ml)

Art. 262023

EAN 4043641 26019 4



Perfume oil vanilla (10 ml)

Art. 262024

EAN 4043641 26020 0



Perfume oil rose (10 ml)

Art. 262025

EAN 4043641 26021 7



Aromaessenz Grapefruit (10 ml)

Art. 262027

EAN 4043641 26023 1



Aromaessenz Limette (10 ml)

Art. 262028

EAN 4043641 26024 8

- 100 % essential oil or composition
- Bottle contents: 10 ml
- Packed in a decorative carton
- Made in Germany

- Trilingual (German, English and French)
- 26 x 26 x 72 mm (LxBxH)
- Suitable for fragrance lamps/aroma diffusers:
e.g. from AL-305, AL-1 300 WS, AL-150 WS, AL-400 CL, AL-240 FL

Special shoes, soles and aids

WE MAKE YOU MOBILE!

HOLISTIC QUALITY SHOE CONCEPT

Special shoes for short-term and long-term care are essential part of the Promed product range. One of the most extensive collections for fields such as therapy, rehabilitation and postoperative care was able to be developed through continuous further development in cooperation with orthopaedists, medical specialists, rehabilitation centres, geriatric nurses and diabetes specialists. Medical studies document the high standard of individual models.

Optimal functionality, high-quality materials and first-class workmanship on the one hand, and innovative ideas and future-oriented technology on the other hand, constitute this holistic quality shoe concept.



Softstep sole concept post-op shoes

The outsole of bandage shoes usually doesn't look very spectacular. In most cases, it is not intended to because stimulation of the healing process is the main focus. This is the very reason why Promed together with specialists found several new solutions which provide relief for people in need and which may stimulate the healing process.



Front without profile

(Relief from Neuropathies and a shuffling gait)



Metashank

A special inlay relieves forefoot metatarsal pain and releases every step

Cylindrical hollow spaces

distributed over the entire sole considerably reduce weight and allow for excellent cushioning properties, especially at the heel

Wider sole

provides for safe walking and standing

Diamond pattern profile

with the "U" format to prevent slipping indoors and outdoors, for a harmonious stride and best damping performance

20 degree beveled heel

provides for reliable step profile and an ideal foot rolling motion

Outsole:

- grindable
- orthopedically adjustable
- made of high quality PU

THE NEWLY DEVELOPED PROMED SOFT STEP SOLE CONCEPT
COMBINES SAFETY AND COMFORT TO A PARTICULARLY HIGH DEGREE.

AIRMESH – THE 'AIR CONDITIONING' UPPER

(upper-meshline)



Interior

AirMesh with EasyDry

The 3-ply material (EasyDry, CosyFit, Mesh) provides specially in summer for well-conditioned feet. From the medical field is derived EasyDry fibre which

removes moisture quickly away from the foot and dries quickly. Particularly suitable for Allergy sufferers.

Style: TheraLight 1, TheraLight 2-S



Interior

AirMesh + Larina-wool

The 3-ply material (pure new wool, CosyFit, Mesh) provides especially in Winter, for toasty warm feet but do not overheat.

Models: GentleWalk-Hi, GentleWalk-Lo




Its network structure and the 3 mm porous CosyFit-special foam lining makes it breathable, so it has 'air conditioning' and provides a perfect feeling of comfort for your feet.

It comes from the medical field and is free of toxic ingredients. Ideal for Allergy sufferers and particularly suitable for diabetics.



COSYFIT – THE BREATHABLE CORE MATERIAL

(insole masterline)



Interior

Models:

- all Masterline, except Munchen
- all Classicline, except Munchen
- Classic 1 and 2



Consists of breathable CosyFit-Material. This material provides for an especially good wear comfort and better pressure relief at the sensitive places.

The multi-layered material, absorbs moisture of the foot to keep dry, and after wearing dries quickly thus provides for a good foot climate.



SYMBOLS



Ideal for orthopedic use



Washable



Polyethylene Suede-look



Sole for super-comfort is widened in the forefoot, to offer trip protection & slip resistance



Washable at 30°



Cosyfit



Soft foam



Micro-Velour



Terry cloth lining



A corn (also clavius from the Latin word for nail) is a callus formation which is caused by chronic pressure on the skin close to bone with central (nail-shaped) spores directed downwards and which is often painful.



Pain in the forefoot may have various causes. Ailments are mostly caused by inflammation due to wrong footwear and/or overload. Pain may also be caused by foot deformities due to accidents.



The different kinds of inflammation often affect joints, tendons, muscle strands or the periosteum itself. Fasciitis plantaris affects the plantar fascia of the sole.



A hammer toe is the deformation of the toe with a flexing contraction of the end joints of the toe.



Callused skin is caused by pressure or friction which causes minimum inflammation in the skin; in order to protect the cellular layers below, the skin is locally stimulated to thicken the outer epidermis.



The diabetic foot syndrome, colloquially also referred to as „diabetic foot“, is a syndrome in connection with diabetes mellitus which mostly occurs in patients with type 2 diabetes.



An arthritic foot is the inflammation protection resulting from inflammation in the joints of the foot or toes.



When suffering from a splayfoot, the bones of the middle foot spread apart.



Deformities of the toe are pathological changes in shape which are either innate or have been developed over the years. Increased mechanical pressure may cause painful calluses and inflammation.



In medical science, Hallux valgus is the x-shaped buckling of the large toe in the base joint. This disease almost always occurs together with a splayfoot.



The multi-layered **COSYFIT** material absorbs the moisture of the foot well & dry quickly after wear, and thus provides for a perfect foot climate.



Particularly suitable for pressure sensitivity of the feet, especially for **diabetics** and all those who are maintaining and healing foot problems and need support.



Particularly suitable for rheumatic symptoms such as **rheumatic** foot changes, claw and hammer toes, hallux rigidus.



The newly developed Promed **SOFTSTEP**-sole-concept combines safety and comfort to an unprecedented level.



The upper material is based on a network structure and the porous **Coys-Fit** breathable foam lining offers a special 'air-conditioning'. Particularly suitable for Allergy sufferers and diabetes.



Promed **'washable-Insoles'** can easily be removed at any time, and washed at 30 degrees. This will prevent unpleasant odours.



Natural new **sheep's wool** ensures warm and cosy feet and wonderfully soft wearing comfort in the whole shoe.



The upper is made from a perforated **stretch micro suede material** for extra breathability.



The shoes are produced in the EC and undergo quality checks.



This ensures utmost quality. Insoles are the substructure of shoes. They form a permanent, inherently stable construction even without the outsole which therefore can be provided with other features.



The resilient material optically resembling suede allows for a healthy breathable climate in your shoe. It is abrasion-resistant, tear-resistant and yet easy to clean.

General information on the range of Promed shoes:



Machine wash
at 30 °C



Not
dryer-proof



Do not
iron



Do not
dry-clean



Do not
use bleach

The shoes can be disinfected using any commercially available or common household disinfectant.

Overview of shoes and indication

Shoe	Gentle Walk-Hi	Gentle Walk-Lo	München 2 LXL	München 3	Wallgau	Grainau	Flex 2
Product image							
Art.	510341	510331	510230	521241 (black) 521242 (brown)	520210	520220	520230
HiMiNr. (Pair)	31.03.03.4058	31.03.03.4059	31.03.03.4059	31.03.03.4151	31.03.03.4152	31.03.03.4149	31.03.03.4150
Page inside the promed catalogue 2023	61	61	62	62	58	58	57

Main fields of application

With general problems in the area of the toes and forefoot				✓	✓	✓	✓✓
Ideal for diabetes prophylaxis	✓✓	✓✓		✓✓	✓✓	✓✓	✓✓
Ideal for rheumatics	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
Ideal with neuropathies	✓	✓	✓	✓	✓	✓	✓
Ideal for wheelchair users	✓✓	✓✓	✓	✓✓	✓✓	✓✓	
With almost all foot shapes, e.g. hallux, hammer toes, osteoarthritis, swollen feet			✓✓		✓	✓	✓
With inflammations, swellings or bandaged feet			✓		✓	✓	✓✓
With foot deformations			✓✓	✓✓	✓✓	✓✓	✓✓
For use after operations					✓		
With rheumatic feet	✓✓	✓✓		✓	✓	✓	✓
For strongly burdened feet (e.g. overweight)							
Wounds in the foot area, e.g. traumatic injuries					✓	✓	✓
With pressure sensitivity, inflammations, swellings and lymphatic feet	✓✓	✓✓			✓		✓
Wounds in the area of the forefoot and on the sole of the foot caused by diabetes mellitus							
Wounds in the area of the forefoot and on the sole of the foot caused by other indications, e.g. traumatic injuries							
Wounds in the rear part of the foot and on the sole of the foot caused by diabetes mellitus							
Wounds in the rear part of the foot and on the sole of the foot caused by other indications, e.g. traumatic injuries							

✓ = suitable ✓✓ = particularly suitable

Theramed D1	Theramed D3	Pantomobil D	TheraLight 1	Sanicabrio DXL	Sanicabrio	Sanicabrio LXL	Sanisoft D	Sanisoft
								
510160	510180	510140	510311	510130	521120	510150	510170	521110
31.03.03.4042	31.03.03.4041	31.03.03.4043	31.03.03.4059	31.03.03.4045	31.03.03.4020	31.03.03.4046	31.03.03.4047	31.03.03.4022
63	63	64	64	59	60	59	60	61

✓✓	✓✓	✓✓	✓✓	✓✓	✓	✓	✓	✓
					✓✓			✓✓✓
✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓	✓✓
✓	✓	✓	✓	✓	✓✓	✓	✓	✓✓
				✓	✓✓	✓	✓	✓✓
✓	✓✓	✓	✓	✓✓		✓✓		
✓✓	✓✓			✓✓		✓✓	✓✓	
✓✓	✓	✓		✓✓		✓✓	✓✓	✓
✓	✓			✓		✓	✓	
					✓✓			✓✓
✓			✓	✓		✓	✓	
✓	✓	✓		✓✓	✓	✓✓	✓	✓✓
✓✓	✓✓	✓✓		✓	✓	✓✓	✓✓	✓
✓✓	✓✓		✓	✓✓		✓✓	✓✓	
				✓✓	✓✓	✓✓	✓✓	✓✓
	✓✓	✓✓		✓		✓	✓	
		✓✓		✓✓			✓	

SHOE DIMENSIONS PROMED MASTERLINE

Models			SaniCabrio DXL, LXL, SaniSoft D			Theramed D 1			Theramed D 3, Pantomobil D			München 2 LXL		
Design			close toe / close heel			open toe / close heel			open toe / open heel			closed shoe		
English Sizes *)	German / French Sizes	Foot length in mm	Outsole in mm	Insole in mm	bale dimension in mm „J width“	Outsole in mm	Insole in mm	bale dimension in mm „J width“	Outsole in mm	Insole in mm	bale dimension in mm „J width“	Outsole in mm	Insole in mm	bale dimension in mm „J width“
3½	36	229	253	238	232	249	239	232	243	239	232	243	234	232
4½	37⅓	237	256	242	237	253	243	237	249	244	237	253	238	237
5	38	241	264	249	240	260	250	240	253	248	240	256	242	240
6	39⅓	250	273	258	245	264	255	245	260	257	245	264	249	245
6½	40	254	277	261	247	269	259	247	264	262	247	273	258	247
7½	41⅓	262	285	268	252	277	268	252	273	271	252	277	261	252
8	42	267	295	276	255	285	275	255	280	278	255	285	268	255
9	43⅓	275	298	280	260	295	283	260	285	281	260	295	276	260
9½	44	279	307	290	262	298	288	262	295	289	262	298	280	262
10½	45⅓	288	312	294	267	307	296	267	298	295	267	307	290	267
11	46	292	315	300	270	312	301	270	307	303	270	312	294	270
12	47⅓	301	324	311	275	318	310	275	312	307	275	315	300	275

SHOE DIMENSIONS PROMED MESHLINE

Models			Gentlewalk Hi SaniCabrio Wallgau, Grainau			Gentlewalk Lo SaniSoft München 3			Theralight I			Flex 2		
Design			close toe, close heel			close toe, close heel			open toe, close heel			close toe, open heel		
English Sizes *)	German / French Sizes	Foot length in mm	Outsole in mm	Insole in mm	bale dimension in mm „G width“	Outsole in mm	Insole in mm	bale dimension in mm „G width“	Outsole in mm	Insole in mm	bale dimension in mm „G width“	Outsole in mm	Insole in mm	bale dimension in mm „J width“
3½	36	229	247	238	222	247	238	222	247	235	222	247	238	222
4½	37⅓	237	254	242	227	254	242	227	250	240	227	254	242	227
5	38	241	263	249	230	259	246	230	259	247	230	263	249	230
6	39⅓	250	268	254	235	273	258	235	263	249	235	268	254	235
6½	40	254	279	261	237	279	261	237	273	258	237	279	261	237
7½	41⅓	262	284	268	242	284	268	242	281	270	242	284	268	242
8	42	267	293	276	245	291	275	245	291	278	245	293	276	245
9	43⅓	275	299	280	250	299	280	250	291	278	250	299	280	250
9½	44	279	306	290	252	306	290	252	291	278	252	306	290	252
10½	45⅓	288	311	294	257	311	294	257	311	300	257	311	294	257
11	46	292	320	303	260	313	300	260	313	302	260	320	303	260
12	47⅓	301	321	311	265	321	311	265	321	312	265	321	311	265

*) 1 englischer Stich = ⅓ inch, also 0,84666667 cm

Art. 540131

promed

Pedibelle Anna

The fancy ankle-high boot for women



**Especially for the lady with pressure-sensitive feet.
Beneficial for joint and back problems.**

- Comfortable, cuddly soft, pleasantly wearable
- 2-Strip-Easy-Fit-System for perfect foot adjustment
- "Vita" anti-slip outsole with continuous strip for secure walking and standing
- Feminine heel for relieving the strain on the Achilles tendon
- Perfect heel fit thanks to incorporated heel cap
- Heel area tapered by 20° for ideal roll-over behaviour
- Sole with cylindrical cavities from the toes into the heel area for best cushioning
- Sizes: 36-42, Comfortable width G



Art. 540112 (black)

Art. 540111 (brown)

promed

Pedibelle Eleine

The neat low shoe for women



**Especially for the lady with pressure-sensitive feet.
Beneficial for joint and back problems.**

- Comfortable, cuddly soft, pleasantly wearable
- Shoe can be opened wide and put on quickly with the 1-Strip-Easy-Fit-System
- "Vita" anti-slip outsole with continuous strip for secure walking
- Feminine heel for relieving the strain on the Achilles tendon
- Perfect heel fit thanks to incorporated heel cap
- Heel area tapered by 20° for ideal roll-over behaviour
- Sole with cylindrical cavities from the toes into the heel area for best cushioning
- Sizes: 36-42, Comfortable width G

promed

Art. 540142

Pedibelle Alexander

The ankle-high boot for men

promed



For all those who are seeking gentle support and yet want great comfort for sensitive feet. Beneficial for strained feet and pressure sores.

- Comfortable and pleasantly wearable
- 2-Strip-Easy-Fit-System for perfect adjustment to virtually any foot shape
- "Vitus" wide anti-slip outsole with wide heel for secure walking
- Good heel fit thanks to incorporated heel cap
- Widened forefoot for comfortable fit and secure footing
- Slight heel height for relieving the strain on the Achilles tendon
- Heel area tapered by 20° for ideal roll-over behaviour
- Sole with cylindrical cavities on the heel for best cushioning
- Sizes: 40-46, Comfortable width G



Art. 540122

Pedibelle Jonas

The comfortable low shoe for men

promed



For all those who are seeking relief for their sensitive feet. Beneficial for pressure sores, swelling and vein problems.

- Comfortable, cuddly soft, pleasantly wearable
- Shoe can be opened wide and put on quickly with the 1-Strip-Easy-fit-System
- "Vitus" anti-slip wide outsole with wide heel for secure walking
- Good heel fit thanks to incorporated heel cap
- Widened forefoot for comfortable fit and secure footing
- Slight heel height for relieving the strain on the Achilles tendon
- Heel area tapered by 20° for ideal roll-over behaviour
- Sole with cylindrical cavities on the heel for best cushioning
- Sizes: 40 - 46, Comfortable width G



Art. 520230 (in pairs) | HiMiNr.: 31.03.03.4150

Flex 2

Unique bandage shoe shaped as slip-on shoe

promed



Heel strap
removable




- Comfortable slip-on shoe for indoor and outdoor use
- Variable heel strap for 3 possibilities to wear
- Improvable fit thanks to variable heel and instep strap
- Functioning of the heel strap: pivotable or removable due to two press studs, heel width adjustable due to hook and loop fastener
- Broadened forefoot for comfortable fit and secure footing
- Instep strap with hook and loop fastener for perfect and stepless adjustment of the instep height
- Ergonomically shaped, orthopaedically adjustable outsole with neuropathic zone in the forefoot
- Washable at 30 °C
- Sizes: 36-47, comfortable width J
- Colour: black

Wallgau

Comfortable low bandage shoe

Therapeutic comfort shoe

High-quality,
double insole
for soft
step comfort



Instep hook and
loop fastener for
extra convenience
when slipping
into and out of
the shoe



- Toe box made of a comfortable stretch material
- Particularly suited for pressure-sensitive feet
- Optically appealing upper material made of elastane (stretch) at the forefoot for particularly good breathability; leather optics at the rear
- Ergonomically shaped, orthopaedically adjustable outsole with neuropathic zone in the forefoot

- Stretch material at the forefoot allows for soft and flexible adjustment for almost all foot shapes (Hallux, hammer toes, arthrosis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes prophylaxis
- Washable at 30°C
- Sizes: 36-47, comfortable width G,
- Colour: black



Art. 520220 (in pairs) | HiMi no.: 31.03.03.4149

Grainau

Therapeutic comfort shoe

High-quality,
double insole for
soft step comfort



Instep hook and loop
fastener for extra
convenience when
slipping into and
out of the shoe



- Toe box made of a comfortable stretch material
- Particularly suited for pressure-sensitive feet
- Optically appealing upper material made of elastane (stretch) at the forefoot for particularly good breathability; leather optics at the rear
- Ergonomically shaped, orthopaedically adjustable outsole with neuropathic zone in the forefoot

- Stretch material at the forefoot allows for soft and flexible adjustment for almost all foot shapes (Hallux, hammer toes, arthrosis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes prophylaxis
- Washable at 30°C
- Sizes: 36-47, comfortable width G,
- Colour: black

Sanicabrio DXL

Over ankle, extra special shoe



- With wide-opening convertible – closure system in the forefoot and heel Velcro for easy access
- Soft and flexible adaptation
- Ideal for swollen feet, hammer toes, use after operations and when foot in plaster or bandaged
- Sizes: 36-47; J width
- Color: navy blue

Wide-opening – convertible – locking system in Forefoot and heel part for easy entry



Sanicabrio LXL

Lower cut, extra-wide special shoe



- With wide-opening convertible – closure system in forefoot area
- Velcro heel provides easy access
- Lower strap to the ankle
- Ideal solution for foot problems, e.g. bandaged feet, heavily loaded, pressure-sensitive feet, Hallux and hammer toes
- Sizes: 36-47; J width
- Color: navy blue

Wide-opening - convertiblelocking system in forefoot and heel part for child's play entry



Sanicabrio

Rehabilitation shoe with forefoot and heel strap



- Individual adjustment facility
- Also suitable as a primary care shoe
- Ideal for extreme pressure sensitivity and severe swelling
- Deep rim in the sole for maximum stability for example with foam inserts
- Terry toweling lining
- Sizes: 36-47; Comfortable width G
- Farbe: navy blue



Wide-opening - convertible-locking system in forefoot and heel with velcro for foolproof entry



Sanisoft D

Lower cut special shoe



- With high volume and adjustable to individual needs in both forefoot and instep
- Climate friendly footwear for diabetics and other feet problems
- Easy to put on and take off and with optimal fitting control
- Especially suitable for wheelchairs
- Sizes: 36-47; J width
- Color: navy blue



Convenient opening to closure by Velcro strips in both the instep and heel



Sanisoft

Rehabilitation shoe in a half-shoe design



- 2 velcro fasteners in the forefoot.
- Heel pressure relief cut-out deep below the ankle
- Deep rim in the sole for maximum stability for example with foam inserts
- Terry toweling lining
- Sizes: 36 - 47; Comfortable width G
- Color: navy blue



Art. 510341 (pair) | HiMiNo.: 31.03.03.4.058 (Hi)

Art. 510331 (pair) | HiMiNo.: 31.03.03.4.059 (Lo)

GentleWalk-Hi / GentleWalk-Lo

Boots with AirMesh upper material and pure newe wool lining



- Breathable AirMesh Upper with Larina-wool lining for perfect foot climate
- 2 wide-opening velcro fasteners in the forefoot for easy entry and firm support
- Pressure-free, perfect fit
- Shock absorbing, soft step sole, for everyday walking in and outdoors

- Suitable for Diabetics and allergy sufferers
- Sizes: 36-47; Comfortable width G
- Color: black

München 2 LXL

Stretch suede-look shoe for therapy and rehabilitation



- Forefoot area enhanced
- 4-way stretch material in the forefoot permit soft and flexible adjustment to serve almost all forms of foot problems (Hallux, hammer toes, arthritis, swollen Feet, etc.)
- The stretchy material allows every movement without painful pressure or wear points
- Foot climate-friendly model
- Sizes: 36-47; J width
- Color: black



The Micro velour upper offers smooth, individual support

Unique interchangeable insole

An enlarged forefoot area and 4-way stretch material in the forefoot provide soft and flexible adaptations to serve virtually all forms of foot problems



brown – Art. 521241 (pair) | HiMiNo.: 31.03.03.4151

black – Art. 521242 (pair) | HiMiNo.: 31.03.03.4151

München 3

Low comfortable bandage shoe made of a material combination of stretch and micro velour



Middle and heel areas in micro velour material for a perfect fit



Velcro system for extra comfort with access and exit

Lightweight – shock absorbing, slip-resistant sole

- Instep hook and loop fastener for extra convenience when slipping into and out of the shoe
- Optically appealing upper material made of a perforated stretch Dacron material for or particularly good breathability
- Ergonomically shaped, orthopaedically adjustable outsole with neuropathic zone in the forefoot
- Particularly suited for pressure-sensitive feet

- Stretch material at the forefoot allows for soft and flexible adjustment for almost all foot shapes (Hallux, hammer toes, arthrosis, swollen feet, etc.)
- Ideal for rheumatic feet and for diabetes prophylaxis
- Washable at 30 °C
- Sizes: 36-47, comfortable width G
- Colour: brown, black

Theramed D1

Multifunctional with 2-hook fabric fasteners and a closed heel



Breathable and lightweight upper

- With two velcro straps at the forefoot and instep
- Open toe, instep strap and closed heel
- Especially recommended for problems toe and forefoot, inflammation, swelling and hallux
- Comfortable in the warmer season breathable and easy care upper
- Sizes: 36-47; J width
- Color: navy blue



complete opening by two velcro fasteners



Art. 510180 (pair) | HiMiNo.: 31.03.03.4041

Theramed D3

Multifunctional 2-hook fabric fasteners with open heel



- Three velcro straps for individual adjustment in the forefoot, instep and heel. Heel has open design
- Recommended for pressure sensitivity in heel, but also with inflammation and swelling in the heel and toes
- Climate-friendly model
- Comfortable in the warmer season with breathable upper
- Sizes: 36-47; J width
- Color: navy blue

Three Velcro fasteners for individual adaptation in forefoot, instep and heel



Pantomobil D

Lightweight, multi-functional, wide-opening special shoe in a slipper design



With wide Velcro forefoot securing for foot care and maintenance and with ability to customise



Breathable Cosy Fit Upper (three layers)

- Recommended for pressure sensitivity, swollen or with inflammation of the forefoot and / or heel.
- Foot climate-friendly model, especially pleasant in the warmer months, thanks to the breathability of the Cosy-fit upper (three layers)
- Includes heel strap
- Sizes: 36-47; Comfortable width J
- Color: navy blue



Art. 510311 (pair) | HiMiNo.: 31.03.03.4047

TheraLight 1

Rehabilitation and comfort shoe with ope toe and closed heel



Slipping on the shoe is super easy thanks to hook and loop fastener system

- Breathable AirMesh Upper with EasyDry-foot for prime foot climate
- Velcro closure that opens wide for easy entry
- Pressure-free, perfect fit
- Shock absorbing soft step sole, for indoors and outdoor use
- Suitable for diabetics and allergy sufferers
- Machine washable at 30°C
- Sizes: 36-47; Comfortable width G
- Color: black

Breathable AirMesh upper material with EasyDry lining



Art. see below (pair)

Replacement flexible foam insole (double sole)



Replacement insole (double sole)
For models Wallgau and Grainau
For shoe sizes: 36-47, comfortable width J
Art. 550236

- Insole made of flexible foam 1: Cotton EVA (3.6mm), black
- Insole made of flexible foam 2: Polyester EVA (4mm), black

Art. see below (pair)

Replacement foam insole



Size: 36-47, J width (pairs)
Color: black
Art. 550230 (model: Sanicabrio/DXL/LXL/DS; Sanisoft/D)
Art. 550231 (model: Theramed D1, Theralight 1)
Art. 550232 (model: Theramed D3, Pantamobil D, Theralight 2-S)

- High-quality insole for soft comfort.
- Terry toweling lining
- 4 mm thick EVA insole

Art. 550235 (pair)

Replacement foam insole



- High-quality insole for soft comfort.
- Terry toweling liner
- 4 mm thick EVA insole
- Suitable for each model of the series Promed
- For shoe sizes: 36-47; Comfortable width G
- Color: White

Art. 550240 (pair)

Replacement fur-feel insole



- Easy-care Insoles for sensitive feet
- High-quality, 100% wool layer
- EVA Insole in 4 mm thickness with soft tread comfort
- Suitable for the models: GentleWalk-Hi & GentleWalk-Lo
- For shoe sizes: 36-47; Comfortable width G
- Color: White



HUBER
Fußtechnologie

by **promed**

IT IS IMPORTANT TO BUILD IN THE RIGHT FOUNDATION


Art. 552055

Hallux Valgus

To prevent further divergence of the big toe



HUBER
Fußtechnologie



Excellent for the following symptoms:

- Prophylaxis in hallux valgus inclination
- Stabbing pain in the ball of the foot
- Dull pain without loading
- Local tenderness

Particularly suitable for:

- People with flat feet and / or large toe deformity
- People with hallux valgus problems

Features:


- Moisture-absorbing, quick drying, skin-friendly
- Free of contaminants
- Microfiber soft material with pleasant comfort

Sizes (1 pair): 35-49

Art. 552050

Heel pain insole

Heel enclosure with u-format pad construction



Excellent for the following symptoms:

- Stabbing intermittent pain
- Dull pain when weight pressure applied
- local tenderness

Particularly suitable for:

- Heel spur, heel spur overload
- People in walking and standing occupations
- Allergy sufferers and diabetes

Features:

- Moisture-absorbing, quick drying, skin-friendly
- Pollutant-free micro-fibre soft material with pleasant comfortable feel


Sizes (1 pair): **XS**/35-37, **S**/38-40, **M**/41-43, **L**/44-46, **XL**/47-49

Art. 552040 Joint relief sole 1 (up to 80 kg)

Art. 552045 Joint relief sole 2 (over 80 kg)

Joint relief (</> 80 kg)

For targeted weight-specific pain relief



Suitable for heavy legs



Suitable for tired legs



sweaty feet & burning feet

Excellent for the following symptoms:

- To relieve pains in the spine and joints
- Complaints with the middle foot bones, joints (Metatarsal)
- Foot pad and ankle burning, joint problems (ankle, knee, hip)

Particularly suitable for:

- People in walking and standing occupations
- Painful, burning, too heavily laden joints, tired pressure-sensitive feet

Features:

- Moisture-absorbing, quick drying, skin-friendly
- Pollutant-free micro-fibre soft material with pleasant comfortable feel

Sizes (1 pair):

XS/35-37, **S**/38-40, **M**/41-43, **L**/44-46, **XL**/47-49




Therapeutic Shoes			Preventive Care	Postoperative Therapy	Conservative Therapy	Trauma Care, Fractures	Sprains, Dislocations	Ligament Insufficiency	Plantar Fasciitis	Diabetic Foot, Neuropathy	Wound Care	Children Orthopaedic Care
Post-op shoe												
MedSurg		Page 70			✓		✓				✓	✓
MedSurg Pediatric		Page 70			✓		✓					
MedSurg Pro®		Page 71			✓		✓				✓	✓
AllRound Shoe®		Page 73			✓	✓	✓				✓	✓
WCS® Light		Page 72	✓	✓	✓						✓	✓
WCS®		Page 72	✓	✓	✓						✓	✓
Commodus Open®		Page 73	✓	✓	✓						✓	✓
Off-loading Shoes												
OrthoWedge		Page 74			✓		✓				✓	✓
OrthoWedge Light		Page 74			✓		✓				✓	✓
Relief Dual®		Page 75			✓		✓				✓	✓
Relief Dual® Plus		Page 75			✓		✓				✓	✓
HeelWedge®		Page 76			✓		✓				✓	✓
Cast Shoes & Boots												
Body Armor® Cast Shoe		Page 76					✓	✓				
SlimLine®		Page 77					✓	✓				
SlimLine® Pediatric		Page 77					✓	✓				✓
Balance shoe												
TwinShoe®		Page 78										



			Preventive Care	Postoperative Therapy	Conservative Therapy	Trauma Care, Fractures	Cast Boots	Sprains, Dislocations	Ligament Insufficiency	Plantar Fasciitis	Diabetic Foot, Neuropathy	Wound Care	Children Orthopaedic Care
Orthoses, Walkers & Bandages													
Lower Leg													
	AirTraveler Walker Page 79				✓		✓			✓		✓	✓
	Relief Insert® Walker 2.0 Page 79				✓	✓	✓		✓	✓		✓	
	Relief Insert® Walker 2.0 Achill Page 80				✓	✓	✓		✓	✓		✓	
	Body Armor® Pro Term Page 80				✓		✓					✓	✓
Ankle													
	Body Armor® Ankle Flex Page 81		✓				✓		✓	✓			
Foot													
	TAS® Page 81				✓								
Positioning Devices													
Positioning Splints													
	NightSplint® Page 82		✓		✓		✓				✓		
	Body Armor® Night Splint Page 83		✓			✓					✓		
Positioning Splint													
	Body Armor® Heel Reliever Page 83		✓									✓	✓
	DARCO Posterior Splint Fix Page 82				✓		✓		✓				
Sole systems and accessories													
	Puzzle Insole® Page 80				✓	✓							
	Relief Contour Insole Page 84				✓	✓						✓	✓
	PegContour® Insole Page 85				✓	✓	✓					✓	✓
	PegContour® Insole Uni Page 86				✓	✓	✓					✓	✓
	PegAssist® Insole Page 87		✓		✓	✓						✓	✓
	Replacement soles for WCS®/WCS® Light Page 87		✓		✓	✓						✓	✓

MedSurg Healing Shoe



Compatible with: replaceable insoles (p. 84 - 87)

**Puzzle Insole®, PegContour® Insole,
PegAssist® Insole, Relief Contour Insole,
PegContour® Insole Uni**

for targeted pressure distribution / pressure relief



INDICATIONS

- Used postoperatively for stability and off-loading, e.g. after foot osteotomy procedures

PURPOSE

- Stability, impact and pressure relief after surgical procedures
- Accommodates even the bulkiest dressings

FEATURES

- Breathable mesh-upper
- Patent-protected sole stiffener
- Replaceable insole

- Can be worn on the right and left

Sizes (s.a. p. 88):

Women: WS/34-36.5, WM/37-38.5, WL/39-41

Men: MS/39-41, MM/41,5-43, ML/43.5-45, MXL/45.5-47

Color: Black

You can order the Toe Cap optional.

MedSurg Pediatric Healing Shoe for children



INDICATIONS

- Used postoperatively for stability and off-loading,
- E.g. after foot osteotomy procedures

PURPOSE

- Stability, impact and pressure relief after surgical procedures
- Accommodates even the bulkiest dressings

FEATURES

- Breathable mesh-upper
- Patent-protected sole stiffener
- Replaceable insole
- Replaceable insole

Sizes (s.a. p. 88): till 32

Color: Black

MedSurg Pro® (High)

Healing Shoe

after surgery – ankle-covering



Compatible with: replaceable insoles (p. 84 - 87)

**Puzzle Insole®, PegContour® Insole,
PegAssist® Insole, Relief Contour Insole,
PegContour® Insole Uni**

for targeted pressure distribution / pressure relief



INDICATIONS

- Used postoperatively for stability and offloading, e.g.
- After foot osteotomy procedures

PURPOSE

- Stability, impact and pressure relief after surgical procedures
- Accommodates insole replacement
- In combination with the PegAssist® creates interim solution for the diabetic foot in need of urgent care

FEATURES

- Enlarged heel counter for improved rear foot stabilization
- Low profile – no height compensation necessary

- Can be worn on the right and left

Sizes (s.a. p. 88):

Women: WS/34-36.5, WM/37-38.5, WL/39-41

Men: MS/39-41, MM/41.5-43, ML/43.5-45,
MXL/45.5-47

Color: Black

You can order the Toe Cap optional.

MedSurg Pro® (Low)

Healing Shoe

after surgery – low



Compatible with: replaceable insoles (p. 84 - 87)

**Puzzle Insole®, PegContour® Insole,
PegAssist® Insole, Relief Contour Insole,
PegContour® Insole Uni**

for targeted pressure distribution / pressure relief

**New light,
breathable
upper**



INDICATIONS

- Used postoperatively for stability and offloading, e.g.
- After foot osteotomy procedures

PURPOSE

- Stability, impact and pressure relief after surgical procedures
- Accommodates insole replacement
- In combination with the PegAssist® creates interim solution for the diabetic foot in need of urgent care

FEATURES

- Enlarged heel counter for improved rear foot stabilization
- Low profile – no height compensation necessary

- Can be worn on the right and left

Sizes (s.a. p. 88):

Women: WS/34-36.5, WM/37-38.5, WL/39-41

Men: MS/39-41, MM/41.5-43, ML/43.5-45,
MXL/45.5-47

Color: Black

You can order the Toe Cap optional.

WCS®

Bandage and wound therapy shoe

**INDICATIONS**

- Treatment of diabetic foot syndrome, open wounds and ulcers
- Designed for pressure redistribution
- Postoperative phase after toe and/or forefoot surgery

PURPOSE

- Insole alterations are possible for selective offloading for individual needs
- The rocker style sole produces an ergonomically correct walking pattern

FEATURES

- Includes customizable insoles
- Easy entry and room for bandages

Sizes (s.a. p. 88):

XS/35 - 36 S/37 - 38, S/M/ 39 - 40,
M/41 - 42; L/43 - 44; XL/45 - 46; XXL/47 - 48

Color: Black

Includes: WCS® Insoles For extended care to complement the progressive healing of wounds. A replacement "four-layer" Insole system can be ordered separately. (p. 87)

You can order the Toe Cap optional.

WCS® Light

Bandage shoe (closed)

**INDICATIONS**

- Treatment for diabetic foot syndrome
- Healing shoe for open wounds and mild ulcerations
- Designed for pressure redistribution
- Used for postoperative phase after toe and/or forefoot surgery

PURPOSE

- Insole alterations are possible for selective offloading for individual needs
- The rocker style sole produces an ergonomically correct walking pattern
- Soft and breathable upper provides maximum comfort and leaves extra space for bandages

FEATURES

- Includes customizable insoles
- Replacement "four-layer" Insole available
- Sold in pairs

Sizes (s.a. p. 88): XS/35 - 36 S/37 - 38, S/M/ 39 - 40,
M/41 - 42; L/43 - 44; XL/45 - 46; XXL/47 - 48

Color: Black

Includes: WCS® Light Insoles. For extended care to complement the progressive healing of wounds. A replacement "four-layer" Insole system can be ordered separately (p.87).

AllRound Shoe®

Closed bandage shoe



Compatible with: replaceable insoles (p. 85)

PegContour® Insole

for targeted pressure distribution



INDICATIONS

- Postoperative care
- Wound care, Arthritis

PURPOSE

- Sole allows for easier mobility without compromising stability
- Sandwich sole construct combines firmness for stability and a cushion sole for softer landings
- Shoe upper is a lightweight, breathable nylon fabric that repels water and soil

FEATURES

- Long ridged internal plate promoting weight distribution
- Low profile sole – no height compensation necessary
- Can be worn on the right or the left

Sizes:

XS/37-38; S/39-40; M/41-42; L/43-44, XL/45-46, XXL/47-48

Color: Black

Commodus Open®

Open bandage shoe



Compatible with: replaceable insoles (p. 84 - 87)

PegAssist® Insole, PegContour® Insole Uni

for targeted pressure distribution / pressure relief



INDICATIONS

- Medial lesions on the MTP I joint, lateral lesions on the MTP V joint and the whole dorsal and plantar forefoot area
- For postoperative treatment and wound management
- Arthritis

PURPOSE

- More stability thanks to high bootleg and adjustable bootleg width
- Heel cap pulled well forward for improved guidance and stability of the heel
- Generous, angular insole cut protects the toes
- Simplified roll-over behaviour thanks to repositioned ball roll
- Wearing very large bandages possible thanks to wide opening

FEATURES

- Individual customisation options, high wearing comfort
- Easy to slip into the shoe as it can be opened completely
- Seamless comfort zone, breathable nylon fabric, workable
- Shock-absorbing EVA outsole, interchangeable insole
- No balance shoe required thanks to the low design height
- The bottom of the shoe is stiffened throughout
- Can be worn on the left or right foot

Sizes (s.a. p. 88): XS/33.5 - 37; S/37.5 - 39.5; M/40 - 41.5; L/42 - 44, XL/44.5 - 47

Color: Black

OrthoWedge Light Off-loading Shoe



Compatible with: replaceable insoles (p. 85 - 86)

PegAssist® Insole, PegContour® Insole Uni
for targeted pressure relief
TwinShoe® for height adjustment (p. 78)



INDICATIONS

- Postoperative after hallux valgus, hammer toe / claw toe surgeries
- Healing of forefoot ulceration

PURPOSE

- Pressure redistributed to the midfoot and rearfoot providing specific relief and protection to the forefoot
- Compression EVA outsole: Abrasion resistant with high cushioning properties

FEATURES

- Plus on comfort: Breathable upper material and removable insole
- Light compression EVA outsole combined shock absorption and stability

- Design of the relief wedge (15°) redistributes pressure to midfoot and rearfoot
- Wide opening and variable fastenings ensure easy access
- Wide Velcro strap with medial and lateral deflection for additional rearfoot support
- Easy to fit on the right or the left

Sizes (s.a. p. 88):

XS/33.5-37; S/37.5-39.5; M/40-41.5;
L/42-44; XL/44.5-47

Color: black

OrthoWedge Off-loading Shoe



Compatible with: replaceable insoles (p. 85 - 86)

PegAssist® Insole, PegContour® Insole Uni
for targeted pressure relief
TwinShoe® for height adjustment (p. 78)



INDICATIONS

- Used postoperatively after procedures to correct Hallux valgus, hammer toes / tailor's bunions
- Forefoot ulcerations

PURPOSE

- Pressure displaced to the midfoot and rear foot
- Especially for the relief and protection of the forefoot

FEATURES

- Replaceable insole also available
- Can be worn on the right and left

Sizes (s.a. p. 88):

XS/33.5-37; S/37.5-39.5; M/40-41.5;
L/42-44; XL/44.5-47

Color: black

You can order the Toe Cap optional.

Relief Dual® Off-loading Shoe



Compatible with: replaceable insoles (p. 84 - 87)

**Puzzle Insole®, PegContour® Insole,
PegAssist® Insole, Relief Contour Insole,
PegContour® Insole Uni**

for targeted pressure distribution / pressure relief



INDICATIONS

- Used postoperatively for stability and off-loading, e.g. after footosteotomy procedures

PURPOSE

- Stability, shock absorption and pressure relief after surgical procedures
- Negative heel due to the use of the modular insole Puzzle Insole®, the cushioned insole Relief Contour Insole and the honeycomb cushioned insole PegContour® Insole
- Accommodates insole replacement (including diabetic-compliant foot bedding)
- In combination with the PegAssist® creates interim solution for the diabetic foot in need of urgent care

FEATURES

- Optimal impact relief – slip-resistant rubber sole – light-weight breathable upper
- Rigid shank throughout entire inner-sole
- Innovative dual-sole technology
- Can be worn on the right and left

Sizes (s.a. p. 88):

Women: WS/34 - 36.5, WM/37-38.5, WL/39-41

Men: MS/39-41, MM/41.5-43, ML/43.5-45, MXL/45.5-47

Color: grey

Relief Dual® Plus Foot-part off-loading shoe with high bootleg



Compatible with: replaceable insoles (p. 84 - 87)

**Puzzle Insole®, PegContour® Insole,
PegAssist® Insole, Relief Contour Insole,
PegContour® Insole Uni**

for targeted pressure distribution / pressure relief



INDICATIONS

- Postoperative for stabilisation and relief e.g. after foot osteotomies
- Postoperative in the event of Haglund's exostosis removals
- Postoperative in the event of fasciitis plantaris
- For wound therapy

PURPOSE

- Optimal support and maximum stability due to high heel cap
- Extreme opening due to the latest re-entry option
- Reduced ankle movement and improved pressure redistribution after foot operations due to a full-length, rigid stiffener in the shoe bottom
- Reduced dorsal extension moments of the distal ankles ensure a higher degree of forefoot relief
- Multiple fastening solution: Extra Velcro fasteners optimise additional support in the ankle joint
- Optimal shock absorption and reduced application of force
- Durable as the bootleg is sewn to the outsole
- Generous, angular insole cut protects the toes
- Stability and optimal rolling
- Negative heel due to the use of the modular insole Puzzle Insole®, the cushioned insole Relief Contour Insole and the honeycomb cushioned insole PegContour® Insole

FEATURES

- Innovative two-sole technology: wear-resistant, light, shock-absorbing, slipproof
- Breathable material
- Flexible, adjustable heel cap
- Caring for bulky dressings
- Space for inserts thanks to a removable insole (also diabetes-adapted foot bed)
- Low design height – no balance shoe required
- Can be worn on the left or right foot

Sizes (s.a. p. 88):

Women: WS/34 - 36.5, WM/37-38.5, WL/39-41

Men: MS/39-41, MM/41.5-43, ML/43.5-45, MXL/45.5-47, MXXL/47.5-50

Color: grey

You can order the Toe Cap optional.

HeelWedge® Off-loading Shoe



Compatible with: replaceable insoles (p. 85 - 86)

PegAssist® Insole, PegContour® Insole Uni
for targeted pressure relief
TwinShoe® for height adjustment (p. 74)



INDICATIONS

- Treatment of injuries or ulcers around the heel (heel ulcers)
- Postoperative care after procedures involving plantar fasciitis
- Heel off-loader after fractures

PURPOSE

- Selective off-loading and pressure relief of the heel

FEATURES

- Replaceable insole
- Integrated toe cap and water-repellent
- Can be worn on the right and left

Sizes (s.a. p. 88): XS/33.5-37; S/37.5-39.5; M/40-41.5; L/42-44; XL/44.5-47

Color: black

Art. 530330

Body Armor® Cast Shoe Cast Boot



INDICATIONS

- Usable as a walking shoe over cast bandages

PURPOSE

- Convenience – cast protection
- Easy entry and easily fitted to the cast
- Ergonomically designed to facilitate a correct walking pattern

FEATURES

- Can be worn on the right and left

Sizes (s.a. p. 88):

XS/24 cm, S/25,5 cm, M/27 cm,
L/29 cm, XL/32 cm, XXL/34,5 cm

Please note that the sizes are specified in cm
(foot length including cast).

Color: black

SlimLine®

Closed bandage shoe



INDICATIONS

- Fractures, cast treatments

PURPOSE

- Effective shock absorbing sole
- Adjustable closure technology for bulky bandages or casts

FEATURES

- Weather-resistant upper material
- Interchangeable insole

- Can be worn on the right and left

Sizes (s.a. p. 88):

XS/33.5-37; S/37.5-39,5; M/40-41,5;
L/42-44; XL/44,5-47

Color: black

Farbe blue: Art. 530430 | HiMiNr.: 31.03.03.4025

Farbe pink: Art. 530433 | HiMiNr.: 31.03.03.4025

SlimLine® Pediatric

Cast and bandage shoe for children



INDICATIONS

- Fractures, cast treatments

PURPOSE

- Ideal protection for the cast
- From wear and tear

FEATURES

- Weather-resistant upper material
- Removable insole
- Can be worn on the right and left

Sizes (s.a. p. 88): S/bis 21; M/22-26; L/27-32

TwinShoe® Balance Shoe

Compatible with

**OrthoWedge, OrthoWedge Light,
HeelWedge®, WCS®,**



INDICATIONS

- To compliment a healing shoe (such as a forefoot off-loading shoe)
- Equalizing leg lengths

PURPOSE

- To ensure a safe and ergonomically correct walking pattern when worn together with healing shoes or off-loading shoes

FEATURES

- Sole height: 36 mm
- Removable insole

- Can be worn on the right and left

Sizes (s.a. p. 88):

XS/35-36 S/37-38, S/M/39-40,
M/41-42; L/43-44; XL/45-46; XXL/47-48

Color: black

You can order the Toe Cap optional.

Toe cap

- 1) for MedSurg, MedSurg Pro®, OrthoWedge, TwinShoe® (black)
- 2) Relief Dual®, Relief Dual® Plus, Relief Insert® Walker 2.0/Achill (grey)
- 3) WCS®



INDICATIONS

- Protects the toes from exposure
- Protects from cold and dirt

PURPOSE

- Protects the toe area from external contaminants
- Increased wearing comfort

FEATURES

- Washable by hand
- Can be worn on the right and left

When ordering, please state your shoe size.

for MedSurg — Art. 530370 + Size:

Women: WS/34-36.5; WM/37-38.5; WL/39-41
Men: MS/39-41; MM/41.5-43; ML/43.5-45; MXL/45.5-47

for MedSurg Pro® — Art. 530395 + Size:

Women: WS/34-36.5; WM/37-38.5; WL/39-41
Men: MS/39-41; MM/41.5-43; ML/43.5-45; MXL/45.5-47

for OrthoWedge — Art. 530360 + Size:

XS/33.5-37; S/37.5-39.5; M/40-41.5; L/42-44; XL/44.5-47

for TwinShoe® — Art. 530390 + Size:

XS/33.5-37; S/37.5-39.5; M/40-41.5; L/42-44; XL/44.5-47

for Relief Dual® — Art. 530396 + Size:

Women: WS/34-36.5; WM/37-38.5; WL/39-41
Men: MS/39-41; MM/41.5-43; ML/43.5-45; MXL/45.5-47

for Relief Insert® Walker 2.0/ Achill

— Art. 532405 + Size:

Women: WS/34-36.5; WM/37-38.5; WL/39-41
Men: MS/39-41; MM/41.5-43; ML/43.5-45; MXL/45.5-47

for WCS® — Art. 530397 + Size:

XS/35-36; S/37-38, S/M/39-40, M/41-42; L/43-44;
XL/45-46; XXL/47-48

AirTraveler Walker Stabilising Orthosis



INDICATIONS

- Acute treatment of sprains of the ankle joint
- Stress fractures
- Fractures of the foot and ankle joint
- Fractures of the metatarsal bones
- Achilles tendon injuries
- Forefoot, midfoot and heel injuries

PURPOSE

- Stabilization & immobilization of the ankle joint
- Prevents rotation, supination and pronation

FEATURES

- Available in a high and low version
- Can be worn on the right and left

Sizes (s.a. p. 88):

XS/bis 36; S/37-38; M/39-41; L/42-44; XL/45-47

Color: Black

Art. 531520 | HiMiNr.: 23.06.01.0070


NEW

DARCO

Relief Insert® Walker 2.0 Lower leg and foot orthosis



Including rucksack
and optional
forefoot strap
for additional
support



INDICATIONS

- TMT 1 arthrodesis (Lapidus arthrodesis)
- Metatarsal fracture / hallux valgus / hallux rigidus / hammer toe and claw toe surgeries
- Fractures of the phalanges
- After rearfoot operations
- In the event of ankle ligament ruptures
- Postoperative after Haglund's exostosis removals (in the event of intact Achilles tendon)
- Preoperative and postoperative / post-traumatic immobilisation of the lower leg/foot region, stable foot or ankle fractures, severe ankle sprains
- Upper ankle arthrodesis

PURPOSE

- For postoperative mobilisation and positioning
- For temporary immobilisation

FEATURES

- Hygienic thanks to the removal of the shoe and interchangeable lining
- Automatic width adjustment
- Soft heel area due to recess in the splint
- Low weight and streamlined design
- Simple wound control thanks to ease of removal and application
- Physiological, safe gait pattern
- Easy to attach and remove, individually adjustable thanks to integrated pump
- Sustainable therapy concept in combination with the foot part relief shoe Relief Dual® Plus
- Can be worn universally on the left or right foot

Sizes (s.a. p. 88):

XS/37-38.5, S/39-41, M/41.5-43, L/43.5-45, XL/45.5-47


Color: blue/grey You can order the Toe Cap optional.

Relief Insert® Walker 2.0 Achill

lower leg and foot orthosis

NEW

DARCO



DAW (DARCO Achilles Wedges)
heel wedges

INDICATIONS

- Achilles tendon ruptures
- Postoperative in the event of Achilles tendon suture
- Postoperative in the event of Haglund's exostosis removals

PURPOSE

- For postoperative mobilisation and positioning with Achilles tendon relief
- For temporary immobilisation

FEATURES

- Hygienic thanks to the removal of the shoe and interchangeable lining
- Automatic width adjustment
- Soft heel area due to recess in the splint
- Low weight and streamlined design

- Simple wound control thanks to ease of removal and application
- Physiological, safe gait pattern
- Easy to attach and remove
- Individually adjustable thanks to integrated pump
- Sustainable therapy concept in combination with the foot part relief shoe Relief Dual® Plus
- Can be worn universally on the left or right foot

Sizes (s.a. p. 88):

XS/37-38.5, S/39-41, M/41.5-43, L/43.5-45, XL/45.5-47

Color: blue/grey You can order the Toe Cap optional.

Body Armor® Pro Term

Feet stump orthosis

DARCO



INDICATIONS

- Diabetic foot syndrome
- PAD (peripheral arterial disease)
- Trauma
- Postoperative for treatment after Chopart amputation or Lisfranc amputation

EFFECT

- Rear loading due to wedge-shaped outsole and inner sole in dorsal extension.
- Stabilisation of the ankle joints by bootleg overlapping the ankle prevents a deformation of the stump
- Blue inlay with integrated pump compresses and thus keeps the stump shape.

FEATURES

- Perfect fit due to air compression – easy to pump up and let out the air via built-in system
- Wedge-shaped sole system for optimal rear loading
- Wide opening makes it easy to put on
- Practical fastening system
- Can be worn on the left and right foot

Sizes: M/16 cm; L/19 cm (Stump length)

Colour: black

TAS® Toe Alignment Splint

Toe Alignment Splint



INDICATIONS

- Postoperative toe alignment after:
- Hallux valgus
- Tailor's bunion
- Hammer toe

PURPOSE

- Strap loops guide toes into proper alignment
- Postoperatively redressing bandage instead of tape bandages

FEATURES

- Soft material is skin friendly
- Can be worn on the right and left

Size: Universal

Adjustable to all foot sizes

Color: white

Replacement straps (white): Item 530351

Body Armor® Ankle Flex

Ankle Bandage



INDICATIONS

- Osteoarthritis of the upper and lower ankle joint
- Rheumatic diseases
- Articular effusion and swelling caused by osteoarthritis or arthritis
- Instability

PURPOSE

- Stabilization and protection of the ankle joint
- Increase in circulation due to bilaterally incorporated silicone pads

FEATURES

- Plantar reinforcement
- Areas of different strength
- Incorporated silicone pads
- Can be worn on the right and left

Sizes/heel circumference (s.a. p. 88):

XS/up to 19.0 cm; S/19.0-21.0 cm; M/21.0-23.0 cm
L/23.0-25.0 cm; XL/25.0-27.0 cm

Color: blue/grey

(90° Plus): Art. 530287 | HiMiNr.: 23.03.01.0024
(110°): Art. 530285 | HiMiNr.: 23.03.01.0024

NEW

DARCO

DARCO Posterior Splint Fix

Foot positioning orthosis



Two Velcro fasteners optimally fix the foot in place in the orthosis



Incl. absorbent fabric in the Plus Version



Plus-Version

Moderate angulation for patients with limited dorsal extension in the upper ankle

- Two oval recesses ensure good air circulation
- Anatomically high-volume heel shape to avoid pressure marks

INDICATIONS

Generally applicable for conservative and postoperative lower leg, ankle and foot injury therapy:

- Postoperative positioning
- Post-traumatic stabilisation / immobilisation
- Immobilisation after an accident
- Short-term stabilisation (ligament ruptures, dislocation, fractures, wounds, strains, sprains, etc.)
- Postoperative positioning

PURPOSE

- Individual positioning of the patient via two different angles of inclination
- Integrated Velcro strips allow easy wrapping with conventional gauze bandages

FEATURES

- 110° model for Achilles tendon pathologies available in 20° pointed foot positions
- Soft inner padding (PE foam)
- Anatomically optimised fit
- Keeps its shape, low weight
- Two flexible, adjustable straps keep the lower leg/foot in a specific position
- Can be worn universally on the left or right foot

Color: white/blue

Sizes **Foot length*** **Calf length***

S	26.5 cm	34.0 cm
M	29.5 cm	37.0 cm
L	32.5 cm	40.0 cm

* of the splint

Art. 530275 | HiMiNr.: 23.03.01.0007

DARCO

NightSplint®

Splinting System



INDICATIONS

- Shortening of the plantar fascia
- Plantar fasciitis
- Tendinitis of the achilles tendon
- Prevention of equines deformities
- Postoperative immobilization

PURPOSE

- Pain relief due to tendon shortenings and inflammation

- In the region of the heel
- Effective foot positioning for therapeutic healing

FEATURES

- Can be worn on the right and left

Sizes:

S/35-36; M/37-39; L/40-43; XL/44-49

Color: black

Body Armor® NightSplint Splinting System



INDICATIONS

- Shortening of the plantar fascia
- Plantar fasciitis
- Inflammation and weakness of the achilles tendon
- Metatarsalgia
- Equinus and ankle contractures

PURPOSE

- Supports proper positioning of the foot at 90° to the leg
- Manual stretching of the plantar fascia (using the windlass mechanism model)

FEATURES

- Easy on – Easy off
- Easy to use hook and loop fastener
- High patient compliance
- Can be worn on the right and left


Sizes: Small (34-36), Standard (37-47)

Color: white/blue

Toe Plate are including in the delivery.

nubbed: Art. 530340 | HiMiNr.: 11.11.05.0040
smooth: Art. 530341 | HiMiNr.: 11.11.05.0041

Body Armor® Heel Reliever Heel suspension (one-size-fits-all)



nubbed
foam



individually adjustable

smooth
foam



INDICATIONS

- For prevention and healing of pressure sores on the heel, with diabetic foot syndrome and/or decubitus
- 14 points or fewer on the Braden Scale
- Diabetes with neurologically sensory failure symptoms
- Foot pulse not palpable
- Paraplegia
- Long-lasting immobility or coma

EFFECT

- Effective pressure relief by means of heel suspension through shifting of pressure from the heel onto the Achilles tendon and the triceps surae on the calf

FEATURES

- Flexible size adjustment through simple cutting-to-size of an additional foam cushion
- Adjustable hook and loop bands and comfortable fit
- Quick and easy to put on
- Can be worn on the left and right foot

Size: Universal
(25 - 54 cm calf circumference, up to shoe size 47)

Colour: grey-blue


INCLUDED IN THE SCOPE OF SUPPLY

- Forefoot cushion, positioning wedge; Washbag
- Additional foam cushion



Puzzle Insole®

Modular insole for pressure distribution



**Make two out of one -
for individual adaptation!**

The Puzzle Insole® enables plantar pressure relief by simply replacing modules. Both EVA foam soles differ only through the different degrees of shore stiffness. The three replaceable modules (small toe module, big toe module and middle module) can be exchanged individually for different indications. The modular insole for planter pressure redistribution has been developed for plantar pressure redistribution, especially after forefoot and metatarsal operations.

INDICATIONS

Little toe module (narrow): tailor's bunion, Jones fracture

Big toe module (wide): osteotomies: Austin, Scarf, Akin, Youngswick; arthrodeses: Lapidus, metatarsophalangeal joint and metatarsophalangeal joint of the big toe

Middle module: Weil and/or Helal osteotomy I-IV; Morton's neuroma; hammer toe, claw toe correction; PIP arthrodeses

PURPOSE

- Used after forefoot and midfoot surgical procedures
- Module insole system for selective plantar pressure relief



FEATURES

- Interchangeable components
- 2.4° negative wedge distal to proximal
- Insole is designed for use with the right or left boot

Suitable for the following products:

MedSurg, MedSurg Pro®, Relief Insert® Walker 2.0, Relief Dual®, Relief Dual® Plus

Sizes: Women: WS /34-36.5; WM /37-38.5; WL/39-41

Men: MS/39-41; MM/41.5-43; ML /43.4-45;

MXL /45.5-47

Color: black-blue (1 set)

Relief Contour Insole

Cushion Insole



INDICATIONS

- For use following hallux valgus and small toe surgeries (austin, akin, scarf, weil, hohmann, etc.)
- Talipes, valgus/flat foot/splayfoot or vasoneurosis
- In the case of angio-neuropathic changes of the foot where mild foot deformities may persist

PURPOSE

- Contours adapt to the anatomy of the foot allowing plantar pressure to be redistributed throughout the entire surface area
- As much as 90 % redistribution of pressure can be achieved

FEATURES

- Insole made of a high resiliency material
- PU material is customizable and can be altered
- Available for right and left boot
- 2.4 ° negative wedge distal to proximal

Suitable for the following products:

MedSurg, MedSurg Pro®, Relief Insert® Walker 2.0, Relief Dual®, Relief Dual® Plus

Sizes:

Women: WS/34-36.5, WM/37-38.5, WL/39-41

Man: MS/39-41, MM/41.5-43, ML/43.5-45, MXL /45.5-47

Color: blue

PegContour® Insole

Off-loading Cushion Insole

for MedSurg, MedSurg Pro®, Relief Dual® ·
Relief Dual® Plus, Relief Insert® Walker 2.0



INDICATIONS

- For use after hallux valgus and small toe surgeries (Austin, Akin, Scarf, Weil, Hohmann, etc.)
- Talipes valgus, flat foot, splayfoot
- Angioneuropathic changes of the foot in combination with mild foot deformities
- **For selective off-loading for:**
 - Plantar lesions in the forefoot or midfoot area
 - Plantar ulcers/ulcer prevention in the forefoot or midfoot area
 - Pressure points in the forefoot or midfoot area

PURPOSE

- Effective insole system for plantar pressure relief in acute care

- Pressure can be reduced by up to 55 %
- Optimal pressure redistribution is achieved through the contours

FEATURES

- Soft lengthwise and crosswise sole covering made of functional polyester
- Customisable PU material
- Customised fit through the removal of single pegs
- Available for the right and left foot

Sizes:

Women: WS/34-36.5, WM/37-38.5, WL/39-41

Men: MS/39-41, MM/41.5-43, ML/43.5-45, MXL/45.5-47

PegContour® Insole

Off-loading Cushion Insole

for AllRound Shoe®



INDICATIONS

- For use after hallux valgus and small toe surgeries (Austin, Akin, Scarf, Weil, Hohmann, etc.)
- Talipes valgus, flat foot, splayfoot
- Angioneuropathic changes of the foot in combination with mild foot deformities
- Heel wounds
- Heel spurs (plantar fasciitis)
- Healing of ulcers in the heel area (rearfoot ulcers)
- Postoperative use and treatment of the heel

PURPOSE

- Effective insole system for plantar pressure relief in acute care
- Pressure can be reduced by up to 55 %

- Optimal pressure redistribution is achieved through the contours

FEATURES

- Soft lengthwise and crosswise sole covering made of functional polyester
- Customisable PU material
- Customised fit through the removal of single pegs
- Available for the right and left foot

Sizes:

XS/37-38, S/39-40, M/41-42,

L/43-44, XL/45-46, XXL/47-48

Color: black

PegContour® Insole Uni, black
Soft-bed insoles with a honeycomb structure
for OrthoWedge, OrthoWedge Light, HeelWedge®,
Commodus Open®, SlimLine®



INDICATIONS

- After bunion and little toe operations (Austin, Akin, Scarf, Weil, Hohmann etc.)
- Skewed foot, flat foot, splayed foot
- Angio-neuropathic foot changes combined with mild foot deformities
- Healing of ulcers in the heel area (rear foot ulcers)
- Heel spur (plantar fasciitis) and heel wounds
- Postoperative applications and treatments on the heel
- **For selective off-loading for:**
 - Plantar lesions
 - Plantar ulcers / ulcer prevention
 - Contact pressure marks

PURPOSE

- Universal-use insole system for all-over pressure redistribution and selective pressure relief
- Pressure reduction of up to 45 % possible
- Optimal pressure redistribution is achieved by means of adapted contours

FEATURES

- Soft, longitudinally and transversely elastic sole cover made of functional polyester
- Abradable PU material (memory foam)
- Individual adjustment with individually removable pegs
- Can be worn on the left and right foot

Sizes: XS/33.5-37, S/37.5-39.5, M/40-41.5, L/42-47, XL/44.5-47

NEW

PegContour® Insole Uni, blue
Soft-bed insoles with a honeycomb structure
for MedSurg, MedSurg Pro®, Relief Dual®,
Relief Dual® Plus, Relief Insert® Walker 2.0



INDICATIONS

- After bunion and little toe operations (Austin, Akin, Scarf, Weil, Hohmann etc.)
- Skewed foot, flat foot, splayed foot
- Angio-neuropathic foot changes combined with mild foot deformities
- Healing of ulcers in the heel area (rear foot ulcers)
- Heel spur (plantar fasciitis) and heel wounds
- Postoperative applications and treatments on the heel
- **For selective off-loading for:**
 - Plantar lesions
 - Plantar ulcers / ulcer prevention
 - Contact pressure marks

PURPOSE

- Universal-use insole system for all-over pressure redistribution and selective pressure relief
- Pressure reduction of up to 45 % possible
- Optimal pressure redistribution is achieved by means of adapted contours

FEATURES

- Soft, longitudinally and transversely elastic sole cover made of functional polyester
- Abradable PU material (memory foam)
- Individual adjustment with individually removable pegs
- Can be worn on the left and right foot

Sizes:

Women: WS/34.0 - 36.5, WM/37.0 - 38.5, WL/39.0 - 41.0
Men: MS/39.0 - 41.0, MM/41.5 - 43.0, ML/43.5 - 45.0, MXL/45.5 - 47.0

PegAssist® Insole

Pressure relief sole or comb insole
for punctual pressure reduction



INDICATIONS

- Selective offloading for:
 - Plantar lesions
 - Plantar ulcerations
 - Pressure points

PURPOSE

- Effective insole system for plantar off-loading in acute treatments
- Removable pegs allow for individual adjustments
- Used postoperatively or in acute treatments
- Enables offloading of up to 50 %

FEATURES

- Insole is designed for use with the right or left boot

BLUE: MedSurg, MedSurg Pro®, Relief Dual®, Relief Dual® Plus, Relief Insert® Walker 2.0

(Sole with blue honeycombs – art. 551010 + size)

Women: WS/34.0 - 36.5, WM/37.0 - 38.5, WL/39.0 - 41.0

Men: MS/39.0 - 41.0, MM/41.5 - 43.0, ML/43.5 - 45.0, MXL/45.5 - 47.0

BLACK: For OrthoWedge, OrthoWedge Light, HeelWedge®, Commodus Open®, SlimLine®

(Sole with black honeycombs – art. 551020 + size)

XS/33.5 - 37.0, S/37.5 - 39.5, M/40.0 - 41.5,

L/42.0 - 44.0, XL/44.5 - 47.0

Color: black or blue

Art. 530211 (set of 4 items) | HiMi no.: for bandage shoes 31.99.99.0999

4-Layer Insole System for WCS®/WCS® Light



4-Layer Insole System

INDICATIONS

- Treatment and management of diabetic foot syndrome
- Treatment for healing of open wounds, ulcers and ulcerations
- Postoperative phase after toe, forefoot and midfoot surgical procedures
- Select customizable pressure relief

PURPOSE

- Insole system for selective wound off-loading
- Progressive care for the diabetic patient

FEATURES


- Precise alteration options for offloading
- Multi-layered insoles with varying degrees of hardness for individualized adjustments
- Available for right and left foot
- Can be easily cut to size: sterile film, marker and disposable scalpel for customisation

Sizes: XS/35-36, S/37-38, S/M/39-40, M/41-42; L/43-44, XL/45-46, XXL/47-48


SIZE CHARTS DARCO

All measurements are defined in cm and correspond to the insole.
DARCO shoe sizes represent a general reference only. The choice of shoe size is significantly influenced by the volume of the bandages and size of cast.


AllRound Shoe®




Commodus Open®, OrthoWedge, Ortho Wedge Light, HeelWedge®, SlimLine®, TwinShoe®




MedSurg, MedSurg Pro®, Relief Dual®, Relief Dual® Plus, Relief Insert® Walker 2.0/Archill



AirTraveler Walker




WCS® Light, WCS®




SIZE CHARTS DARCO

All measurements are defined in cm and correspond to the insole.
DARCO shoe sizes represent a general reference only. The choice of shoe size is significantly influenced by the volume of the bandages and size of cast.


Body Armor® Cast Shoe




Body Amor® Sleeve



Body Amor® Ankle Flex




Body Armor® Heel Reliever




SIZE CHARTS DARCO CHILDREN




OrthoWedge Pediatric



MedSurg Pediatric



SlimLine® Pediatric



All measurements are given in cm and correspond to the size of the insole.

**Art. 198001****Modellage Set Ceramic****2**

Professional ceramic modelling set for removing, reworking, modelling UV gel and acrylic, beautifying, polishing. Quick and effective working; durable material, low heat development



Contents Art. 198604	Ceramic bit, cylinder, rough, green
Art. 198603	Ceramic bit, large cone, fine, red
Art. 198606	Ceramic milling cutter, long cone, medium, blue
Art. 198607	Ceramic milling cutter, pointed, fine, red
Art. 198133	Silicone polisher, pointed, fine
	Including bit box, high size

Art. 198002**Modellage Set Blue mix****2**

Professional combination set consisting of ceramic and blue mix abrasive wheels. For matting, shortening, removing, reworking, polishing gel and acrylic; durable material, low heat development



Contents Art. 198480	Hard alloy bit, cylinder, blue mix series
Art. 198481	Hard alloy bit, large cone, blue mix series
Art. 198606	Ceramic milling cutter, long cone, medium, blue
Art. 198607	Ceramic milling cutter, pointed, fine, red
Art. 198133	Silicone polisher, pointed, fine
	Including bit box, high size

Art. 198003**Modellage Set Rainbow****2**

Professional hard alloy bit set with tungsten nitride coating. For removing, modelling, filing, finishing work at the artificial nail, polishing; prevents the development of heat.



Contents Art. 198491	Hard alloy bit, large round cone, rainbow series
Art. 198490	Hard alloy bit, cylinder, rainbow series
Art. 198492	Hard alloy bit, narrow cone, rainbow series
Art. 198493	Hard alloy bit, narrow round cone, fine,
Art. 198133	Silicone polisher, pointed, fine
	Including bit box, high size



Art. 198008

Modellage Set HM-Bitset Pro

Professional high-quality nail modellage set – the „must have“ for every nail designer. Removing, shaping, enhancing and quick modelling is no longer a problem.

Content	Art. 198451	HM-Bit, cylinder, comb-toothing (fine/coarse)
	Art. 198425	HM-Bit, narrow cone, violet series
	Art. 198442	HM-Bit, large tipped cone, red series
	Art. 198557	Carbide-Bit, cylinder, cross-toothing
	Art. 198556	Titan-Bit, cylinder, rounded (gold head)

2



Art. 198006

Modellage Set Master

The complete professional kit for every type of nail modellage.

Content	Art. 198410	Hard metal milling cutter, blue series
	Art. 198421	Hard metal milling cutter, yellow
	Art. 198450	Hard metal milling cutter, blue series
	Art. 198540	Diamond milling cutter
	Art. 198031	Sapphire bit, round
	Art. 198551/L	Carbide titanium milling cutter, long
	Art. 198030	Sapphire milling cutter, pointed
	Art. 961910	Design hard box, transparent

2



Art. 198005

Modellage Set Basic

The perfect basic kit for complete nail modellage.

Content	Art. 198421	Hard metal milling cutter, yellow
	Art. 198450	Hard metal milling cutter, blue series
	Art. 198540	Diamond milling cutter
	Art. 198551	Carbide titanium milling cutter, short
	Art. 961910	Design hard box, transparent

2



Art. 198014

Personal Care Set

For the complete manicure of natural nails: From the trimming of nails, removal of cuticles, thinning out thickened areas on the nails, removal of irregularities on the body of the nail.

Content	Art. 190011	Abrasive disc, sapphire
	Art. 190067	Callus sander, sapphire
	Art. 198540	Diamond cone, small
	Art. 198535	Stone sander, standard
	Art. 198541	Diamond bit, cylinder
	Art. 198031	Sapphire bit, round
	Art. 198030	Sapphire bit, pointed
	Art. 190041	Polishing cone, felt

1

1 = for treat of natural nails

2 = for treat of artificial nails

3 = or treat of the skin

**Art. 198009****1**

Natural Set

For the complete manicure of natural nails: From the trimming of nails, removal of cuticles, thinning out thickened areas on the nails, removal of irregularities on the body of the nail, ending perfectly with a mirror finish.

Content	Art. 190011	Abrasive disc, sapphire
	Art. 198020	Diamond bit, cone
	Art. 198540	Diamond cone, small
	Art. 198535	Stone sander, standard
	Art. 198130	Silicone polisher, fine, yellow
	Art. 198131	Silicone polisher, medium, black
	Art. 198132	Silicone polisher, coarse, green
	Art. 190041	Polishing cone, felt
	Art. 198571	Soft polisher, textile

**Art. 198010****3**

Calluses Set

To remove strong calluses and thickened skin on the hands, feet and other areas of the body such as the elbow, and for the removal of cuticles.

Content	Art. 190021	Abrasive cone, sapphire
	Art. 190067	Callus sander, sapphire
	Art. 198540	Diamond cone, small

**Art. 198016****1****3**

Comfort Set

For all problem zones in the skin, cuticle and nail area. Removes hard skin, calluses, excess cuticles, trims nails, produces a matt finish, smoothes, cleans and polishes.

Content	Art. 198535	Stone sander, standard
	Art. 198530	Stone sander, cylinder
	Art. 198564	Abrasive mandrel, plastic
	Art. 198031	Sapphire bit, round
	Art. 198543	Diamond milling cutter, small, pointed
	Art. 198034	Sapphire bit, long
	Art. 198540	Diamond milling cutter
	Art. 198541	Diamond bit, cylinder
	Art. 190086	Goat hair brush
	Art. 190021	Abrasive cone, sapphire
	Art. 190041	Felt cone buffer
	Art. 190011	Abrasive disc, sapphire
	Art. 198562	Abrasive belt, medium (4 pcs.)
	Art. 961910	Design hard box, transparent

*autoclave (bit sterilisation by a temperature of 121-134 °C and a pressure of ca. 2 bar)

1 = for treat of natural nails **2** = for treat of artificial nails **3** = or treat of the skin

Note: All promed abrasive products have a standard shaft dental diameter of 2.332 to 2.350 mm.

Ceramic Bits

More lightweight than hard alloy bits, durable, efficient

Suitable for removing gel and acrylic, filing, engraving, smoothing, polishing and reshaping as well as removing dead skin. Extremely good treatment of the nails possible, smooth surfaces, pleasant feeling on the nails, quick, easy to use, good running properties due to reduced vibration, very resistant to corrosion, low heat development



Art. 198600

Ceramic milling cutter, large cone, fine, red *

2

Suitable for reworking very thickened nails and callus as well as artificial nails, smooth running properties, extremely hard, durable, sharp. Speed range: up to 30,000 rpm



Art. 198601

Ceramic milling cutter, large cone, medium, blue *

2

3

Suitable for treating very hard, thickened nails and callus as well as artificial nails. Ideal for removing UV gel and acrylic effectively. Also suitable for reworking. Lower weight as compared to hard alloy bits, smooth running properties, extremely hard, durable, sharp. Speed range: up to 30,000 rpm



Art. 198602

Ceramic milling cutter, large cone, rough, green *

2

3

Suitable for treating very hard, thickened nails and callus as well as artificial nails. Ideal for removing UV gel and acrylic effectively. Lower weight as compared to hard alloy bits, smooth running properties, extremely durable, sharp. Speed range: up to 30,000 rpm



Art. 198604

Ceramic milling cutter, cylinder, rough, green *

2

For matting the tips, shortening artificial nails, removing gel and acrylic. Very good cutting ability, pleasant feeling on the nails. Speed range: up to 20,000 rpm



Art. 198603

Ceramic milling cutter, large cone, fine, red *

2

An absolute "all-rounder" in its field. The unique cone shape allows you to remove gel and acrylic (also 3-phase gels) with little risk. Moreover, it provides an optimum basis for carrying out additional work at the modelled nails. Very good cutting ability. Speed range: up to 25,000 rpm



Art. 198605

Ceramic milling cutter, long cone, fine, red *

2

For reworking, softly polishing and reshaping your nails. Speed range: up to 20,000 rpm



Art. 198606

Ceramic milling cutter, long cone, medium, blue *

2

For removing, for easy surface treatment, also possible for shortening. irregular nails. Speed range: up to 20,000 rpm



Art. 198607

Ceramic milling cutter, pointed, fine, red *

2

Narrow bit for work at the nail fold/edge, easy trimming or work below the nails. Speed range: up to 35,000 rpm



Art. 198608

Ceramic milling cutter, pointed, medium, blue *

2

Narrow bit for work at small and narrow spots, within easy reach without scratching the nails. For reworking the modelled nails, for removing residues below the nails. Speed range: up to 35,000 rpm

Blue mix bits

Very durable material, prevents the development of heat

Blue mix series: Hard alloy bits with aluminium nitride coating.
Particularly rapid and speedy treatment of artificial nails.



Art. 198480

Hard alloy bit (with aluminium nitride coating), cylinder, blue mix series *

For matting the tips, shortening artificial nails, removing gel and acrylic. Coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2



Art. 198481

Hard alloy bit (with aluminium nitride coating), large cone, blue mix series *

An absolute "all-rounder" in its field. The unique cone shape allows you to remove gel and acrylic (also 3-phase gels) free of risk. Moreover, it provides an optimum basis for carrying out additional work at the modelled nails. Coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2

Rainbow bits

Very durable material, prevents the development of heat

Rainbow series: Hard alloy bits with tungsten nitride coating
Particularly rapid and speedy treatment of artificial nails.



Art. 198491

Hard alloy bit (with tungsten nitride coating), large round cone *

An absolute "all-rounder" in its field. The unique cone shape allows you to remove gel and acrylic (also 3-phase gels) free of risk. Moreover, it provides an optimum basis for carrying out additional work at the modelled nails. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2



Art. 198490

Hard alloy bit (with tungsten nitride coating), cylinder *

For matting the tips, shortening artificial nails, removing gel and acrylic. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 25,000 rpm

2



Art. 198492

Hard alloy bit (with tungsten nitride coating), narrow cone, fine *

Gets acrylic and (building) gel into shape. Due to its size, it is particularly suited for finishing work. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 30,000 rpm

2



Art. 198493

Hard alloy bit (with tungsten nitride coating), narrow small round cone, fine *

Gets acrylic and (building) gel into shape and files it off very well. Due to its shape, it fits the nail shape very well and allows higher performance. Tungsten nitride coating prevents the development of heat, very durable material. Speed range: up to 30,000 rpm

2

* autoclave (bit sterilisation by a temperature of 121-134 °C and a pressure of ca. 2 bar)

1 = for treat of natural nails 2 = for treat of artificial nails 3 = or treat of the skin

Note: All promoted abrasive products have a standard shaft dental diameter of 2.332 to 2.350 mm.

Titanium Bits / Full Carbide-bits

for extrem demands and long life

These high quality bits from titanium-carbit-alloy heats up less fast and are more pleasant therefore in application for your customers.



Art. 198551

Carbide titanium milling cutter, barrel, short shaft *

Dulling the tips, shortening of artificial nails, for removal of gel and acrylic, extrem durable. Rotational speed: up to 25,000 rpm.

2




Art. 198551/L

Carbide titanium milling cutter, barrel, long shaft *

Dulling the tips, shortening of artificial nails, for removal of gel and acrylic, extrem durable. Rotational speed: up to 25,000 rpm.

2




Art. 198556

Carbide titanium milling cutter, rounded barrel *

Tip blending to the natural nail, dulling the tips, shortening of artificial nails, you can file the nails thinner on the hem at the nail, for removal of gel and acrylic, extrem durable. Rotational speed: up to 25,000 rpm.

2



Art. 198557

Full Carbide-Bit barrel, cross cut *

Special grinder for both running directions, fast removal of the modeling (gel and acrylic), special coating avoid heat development, very fine structure on the nails image, extrem durable. Rotational speed: up to 25,000 rpm.

2

Hard alloy grinders

very durable material

Carbide bits are especially well suited to allowing quick and effective work on artificial nails, and should not be used on natural nails. Typical applications: tip blending, filing down surfaces, filing the underside of the nails, shortening, nail gems.

Hard alloy grinders, lilac series, Spiral-cut - fine



Art. 198455

Hard alloy bit, drops *

Nail and skin smoothing as well as nail art work.

Rotational speed: optimal at 12,000 rpm.; max. up to 20,000 rpm.

2



Art. 198435

Hard alloy bit, long cone *

Facilitates a multi-purpose application area: Cleaning underside of the nails on nail smile line. Leftovers of gel or acrylic are removed and finishing touch is done.

Rotational speed: up to 40,000 rpm.

2



Art. 198425

Hard alloy bit, narrow cone *

Is particularly efficient by s.g. Lifting. Finishing touch and quick removal of gel or acrylic is possible, even in the nail fold. Rotational speed: up to 30,000 rpm.

2



Art. 198445

Hard alloy bit, large cone *

Is an absolute ALL-Rounder in its field. The unique cone-shape offers a risk-free removal of gel and acrylic (even 3-phase gel) without generating heat. It also leaves behind an optimal base to work further on the nail modeling. An extremely durable material. Rotational speed: up to 20,000 rpm.

2



Hard alloy grinders, blue series, Cross cut - rough

Art. 198410

Hard alloy grinder, round *

Particularly suitable to plunge stones in the artificial nail. Rotational speed: up to 40,000 rpm.

2



Art. 198420

Hard alloy grinder, cone *

It brings acrylic and (get-up) gel in shape and files very nicely. Its shape fits wonderfully to the shape of the nail and enables a higher workmanship. No heat generation.

Rotational speed: up to 20,000 rpm.

2



Art. 198430

Hard alloy grinder, small cone *

Bring acrylic and (get-up) gel in shape. Its size is particularly ideal for finishing.

Rotational speed: up to 30,000 rpm.

2



Art. 198450

Hard alloy grinder, barrel *

Very quick and rough removal of finger nails (acrylic and gel). Very low heat generation.

Rotational speed: up to 20,000 rpm.

2

especially well suited
for 3-phase gels

Hard alloy grinders, red series, Cross cut - medium




Art. 198422

Hard alloy grinder, cone *

Optimum removal of any kind of gels - especially suitable for a refill. The conical shape is suitable for finishing the nail modeling for e.g. to bring into the shape.

Rotational speed: up to 30,000 rpm.

2



Art. 198442

Hard alloy grinder, bigfootball *

Thanks to its shape, the small pointed cone is ideal for precise ablation of gel and acrylic also to the delicate points just before the nail fold. Rotational speed: up to 25,000 rpm.

2



Art. 198443

Hard alloy grinder, small football *

Thanks to its shape, the small pointed cone is ideal for precise ablation of gel and acrylic also to the delicate points just before the nail fold. Rotational speed: up to 30,000 rpm.

2



Art. 198452

Hard alloy grinder, rounded barrel *

The rounded cylinder shape provides a safe and superficial removal of gel and acrylic on the entire nail plate without generating heat. Rotational speed: up to 20,000 rpm.

2

Only for treating
fine details

Hard alloy grinders, yellow series, Cross cut - extra fine




Art. 198434

Hard alloy grinder, long *

Its fine surface texture (cross-cut) cleans perfectly on the lower side of the nail. Removal of leftovers of acrylic or gel, final touch of the work. Rotational speed: up to 40,000 rpm.

2



Art. 198421

Hard alloy grinder, cone *

Reworking on nail modeling and removing the leftovers. For perfect cleaning of the underside of the nails. Rotational speed: up to 30,000 rpm.

2



Art. 198461

Hard alloy grinder, flame *

HM-Bit Flame reaches even to the most delicate point due to its shape. To remove adhesive or stuff under the nails made of acrylic and gel. It can also be used to clean the natural nails after re-growth. Rotational speed: up to 40,000 rpm.

2

HM-Bits, Specials



Art. 198451

Hard alloy grinder, barrel, (fine, rough)*

Combination cut (with fine and rough elements); very durable; fast removal without heat development; very smooth surface structure on the nails. Rotational speed: up to 20,000 rpm.

2



Art. 198453

HM bit, narrow cone, variable toothing (fine/coarse), orange series

for finishing the shaping and removing residue, finishing work and quick removal without heat development of gel or acrylic, very smooth surface structure on the nail, very durable material, can also be used in podiatry for removing calluses. Speed range: up to 20,000 rpm.

2



Art. 198423

Hard alloy grinder, cone, middle *

especially for rotation direction left; fast removal of the modeling (gel and acrylic) very smooth surface structure on the nails; very durable. Rotational speed: up to 20,000 rpm.

2



Art. 198470

Special milling cutter with cross cut, Green series, rough

Treatment of very thick dystrophic nails and hard skin.

Rotational speed: optimal at 12,000 rpm.; max. up to 20,000 rpm.

2



Art. 198471

Special milling cutter with cross cut, Orange series, fine

Removing hardened skin and thick calluses as well as nail art work.

With external cooling grooves for cool working in dry and wet Technique.

Rotational speed: optimal at 10,000 rpm.; max. up to 15,000 rpm.

2

Steel cutter – stainless



Art. 197005

Round cutter with cross cut (2,3 mm) *

For removing corns and calluses from the skin. High-efficiency toothing with cross cut.

Rotational speed: optimal at 15,000 rpm.; max. up to 30,000 rpm.

3



Art. 197015

Fissure cutter (1,2 mm) *

For treatment of the nail fold. With cross-cut toothing.

Rotational speed: optimal at 15,000 rpm.; max. up to 50,000 rpm.

1 3



Art. 197021

Hollow cutter (trepan, 2,3 mm) *

For removal of cracked hard skin and calluses. With toothed blade.

Rotational speed: optimal at 4,000 rpm.; max. up to 6,000 rpm.

3



Art. 197025

Rotary tissue punch, grounded (2,7 mm) *

For a gentle removal of deep-set corns / calvi. With ground blade.

Rotational speed: optimal at 4,000 rpm.; max. up to 6,000 rpm.

3



Art. 197030

Nail-Cleaner, short *

For a gentle removal of deposits from the nail fold and cuticles. Hexagonal working part without sharp blades. Rotational speed: optimal at 5,000 rpm.; max. up to 30,000 rpm.

1 2 3



Art. 197031

Nail-Cleaner, long *

For a gentle removal of deposits from the nail fold and cuticles. Hexagonal working part without sharp blades. Rotational speed: optimal at 5,000 rpm.; max. up to 30,000 rpm.

1 2 3

Diamond bits

The fine diamond particle coating allows hygienic and speedy skin and nail work. Great durability when used intensively. Smooths, cleans and removes sticky residues, even in the nail fold. Safe operation in problem areas.



Art. 198505

Nail cutter (diamond coated on Top only)*

1 2

Ideal for shortening nails. The adjacent skin remains untouched. Safety bevel at the edge of the front end. Rotational speed: optimal at 15,000 rpm.; max. up to 20,000 rpm.



Art. 198502

Hard skin abrader, cone (super coarse grit) *

3

Removal of hard callus. With external cooling grooves for cool working in dry and wet Technique. Rotational speed: optimal at 10,000 rpm.; max. up to 15,000 rpm.



Art. 198501

Hard skin abrader, conical (super coarse grit) *

3

Removal of hard callus. With external cooling grooves for cool working in dry and wet Technique. Rotational speed: optimal at 10,000 rpm.; max. up to 15,000 rpm.



Art. 198585

Diamond milling cutter, rounded barrel, big, rough, hard skin abrader* (especially for wet technique)

3

For the professional pedicure. No heat development. For the treatment of skin irregularities and calluses. Can be used on moist and dry skin. Rotational speed: up to 20,000 rpm.



Art. 198586

Diamond milling cutter, rounded barrel, small, rough, hard skin abrader* (especially for wet technique)

3

For the professional pedicure. No heat development. For the treatment of very thick calluses. Can be used on moist and dry skin. Rotational speed: up to 20,000 rpm.




Art. 198020

Abrasive cone, diamond*

1 3

To grind away uneven, grooved nail surfaces. Rotational speed: up to 20,000 rpm.



Art. 198540

Diamond milling cutter, small cone*

1 3

For the removal of cuticles, mat the transitions and the natural nail, to treat difficult edges. Rotational speed: up to 20,000 rpm.



Art. 198541

Diamond milling cutter, barrel*

1 3

For the treatment of thickened nails and the hard edges of corns. Rotational speed: up to 20,000 rpm.



Art. 198542

Diamond milling cutter, football*

1 3

For the grinding away of fungal nails, fungal traces can be removed well with this head. For the removal of cuticles. Rotational speed: up to 20,000 rpm.



Art. 198543

Diamond milling cutter, small football*

1 3

For the grinding away of fungal nails, fungal traces can be removed well with this head. For the removal of cuticles. Rotational speed: up to 20,000 rpm.

*autoclave (bit sterilisation by a temperature of 121-134 °C and a pressure of ca. 2 bar)

1 = for treat of natural nails 2 = for treat of artificial nails 3 = or treat of the skin

Note: All promed abrasive products have a standard shaft dental diameter of 2.332 to 2.350 mm.

Saphir-bits

Sapphire is „softer“ than diamond and ideal for treating and caring for natural nails, thickened nails, hard skin and other foot and nail problems.



Art. 198030

Sapphire milling cutter, pointed

To treat difficult corners and notches and to remove cuticles carefully.

Rotational speed: up to 25,000 rpm.

1 3



Art. 198031

Sapphire milling cutter, round

Use for clearing nail fold as well as removal of callosities To grind away ingrown nails.

Rotational speed: up to 25,000 rpm.

1 3



Art. 198034

Sapphire milling cutter, long

To treat the underside of nails, e.g., to remove dirt or the spores from nail fungi.

Rotational speed: up to 25,000 rpm.

1 3




Art. 190011

Abrasive disc, sapphire

To trim and shape fingernails and toenails. Rotational speed: up to 10,000 rpm.

1



Art. 190027 fine

Art. 190029 middle

Art. 190031 rough

Abrasive cone, sapphire small*

For removal of thickened areas or calluses on the sole of the foot. Can be even used for removal of corns on the toes. Rotational speed: up to 20,000 rpm.

1 3



Art. 190024 fine

Art. 190026 middle

Art. 190028 rough


Art. 190021 rough (not autoclave)

Abrasive cone, sapphire*

To grind away calluses, thickened skin and thickened areas on the nails.

Rotational speed: up to 20,000 rpm.

1 3



Art. 190066 fine

Art. 190064 middle

Art. 190063 rough

Art. 190067 rough (not autoclave)


Hard skin abrader, sapphire*

Especially suitable to remove for stro ng and thick calluses, e.g. on the heels and the balls of the feet. Rotational speed: up to 20,000 rpm.

3

Stonegrinders

The stone abraders are perfect suitable for thinning out heavy thick nails.



Art. 198530

Stone abrader, cylindre

To smooth out unattractive nail bodies, for thickened nails.

Rotational speed: up to 12,000 rpm.

1



Art. 198535

Stone abrader, standard

To smooth out and make the body of the nail matte. As preparation for the application of fake nails.

Rotational speed: up to 12,000 rpm.

1

Sandpapergrinders and mandrels

Sanding cylinders are available with various grain sizes. The cylinders are placed on a mandrel, replacing the good old hand file with a well-rounded modern alternative.



Art. 198564

Abrasive mandrel, plastics, autoclave*

Suitable for the placement of abrasive bands to work on the modeling of artificial nails. Rotational speed: up to 25,000 rpm.

2



Art. 198565/L

Abrasive mandrel, plastics, long, autoclave*

(for 4030/S/SX/SX2, 4040, 5030S, 5040-SX2)*

Suitable for the placement of abrasive bands to work on the modeling of artificial nails. Rotational speed: up to 25,000 rpm.

2



Art. 198561 fine

Abrasive band, fine, 100 pcs.

Suitable for grinding down the nail surface and for finely work.

2



Art. 198562 medium

Abrasive band, medium, 100 pcs.

Especially suitable to removal of gel.

2



Art. 198563 rough

Abrasive band, rough, 100 pcs.

Suitable to removal of acrylic.

2



Art. 198590 fine

Abrasive band ZEBRA, fine, 100 pcs

With „Anti-Dust“ coating. Suitable for grinding down the nail surface and for finely work.

2



Art. 198591 medium

Abrasive band ZEBRA, medium, 100 pcs.

With „Anti-Dust“ coating. Especially suitable to removal of gel.

2



Art. 198592 rough

Abrasive band ZEBRA, rough, 100 pcs.

With „Anti-Dust“ coating. Suitable to removal of acrylic.

2



Art. 198580

Abrasive mandrel, rubber (for abrasive cap 13 mm)

Abrasive mandrel for placement of sand paper-caps for the processing of calluses. Rotational speed: up to 15,000 rpm.

3



Art. 198588

Abrasive cap (13 mm), medium, 100 pcs.

Professional abrasive cap for the pedicure. For the removal of heavy and thick calluses, callosities on ball and heel.

3



Art. 198581

Abrasive mandrel, rubber (for abrasive cap 10 mm)

Abrasive mandrel for placement of sand paper-caps for the processing of calluses. Rotational speed: up to 15,000 rpm.

3



Art. 198587


Abrasive cap (10 mm), medium, 100 pcs.

Professional abrasive cap for the pedicure. For the removal of heavy and thick calluses, callosities on ball and heel.

3

Polishing Bits

Polishing, removing shine, glossy finishes, matt finishes, smoothing out unevenness in natural nails and much more. The various polishing bits, from silicone buffers and felt bits to soft textile polishers, offer a wide spectrum of nail care applications.



Art. 190041

Felt cone buffer

Smooths the edges and surfaces of the nails for the polished finish to treatment. Rotational speed: up to 10,000 rpm.

1 2




Art. 198133

Silicone polisher, sharp, fine

Use for polishing, creating a shiny surface and matting.

1 2



Art. 198130

Silicone polisher, fine, yellow

Use for polishing, creating a shiny surface and matting.

1 2




Art. 198131

Silicone polisher, medium, black

Use for polishing, creating a shiny surface and matting.

1 2



Art. 198132

Silicone polisher, rough, green

Use for polishing, creating a shiny surface and matting.

Removal of slight irregularities on natural nails. **Buffer compensation!**

1 2

General:

Three different grains sizes for the fine polishing of natural and artificial nails.

Rotational speed: up to 8,000 rpm

(Attention: over this the nails will be discolored with the color of the silicon polisher)



Art. 198572

Soft polisher, textile, incl. polish

Textile incl. polish creates a unique glossy finish on natural nails, as the polishing paste has already been

integrated into this polishing bit. Speed range: up to 10,000 rpm

1 2



Art. 198570


Soft polisher, leather

For polishing of natural nails. In combination with a polishing paste also for artificial nails.

The leather produced by its surface structure a wonderful gloss on the natural nail plate.

Rotational speed: up to 10,000 rpm.

1 2



Art. 198571

Soft polisher, textile

For the glossy polishing of fake and natural nails in combination with a polish paste.

Rotational speed: up to 10,000 rpm.

1 2



Art. 190086

Polishing brush, goat hair

To clean the nails and the materials. You can remove the dust of the nails on hard reach places, like around the cuticles. Brush is made of soft goat hair. Rotational speed: up to 10,000 rpm.

2

*autoclave (bit sterilisation by a temperature of 121-134 °C and a pressure of ca. 2 bar)

1 = for treat of natural nails

2 = for treat of artificial nails

3 = or treat of the skin

Note: All promed abrasive products have a standard shaft dental diameter of 2.332 to 2.350 mm.

Accessories expert



Art. 198915

Bitholder, plastic

space for 13 grinding bits (delivery without bits), white, 100 x 30 mm



Art. 198900

Acrylic-Bitholder

space for 10 grinding bits (delivery without bits), transparent, 95 x 50 mm



Art. 198910

Cleaning brush for bits, brass


for cleaning of grinding bits.



Art. 198920

Cleaning brush for bits, nylon

for cleaning of grinding bits. Especially gentle for the blades of the grinding bits.



Art. 198930

Cleaning brush for bits, nylon

for cleaning of grinding bits. Especially gentle for the blades of the grinding bits.



Art. 940602 (10 pcs.) / Art. 940601 (1 piece)

Promed towels

white, small, 30 x 50 cm, delivery contents: 10 pcs. or 1 piece



Art. 952010

Handbook for professional nailart

This compact, richly illustrated book explains the most important techniques when working with electric files and covers all the most relevant subjects.

68 pages, 80 illustrations, ISBN 3-9809164-0-5



Art. 190051 - white

Chuck - white (for Pedisenso Duo)

fixed chuck

Art. 190052 - yellow

Chucks - yellow (for Pedisenso Duo)

with safety slip clutch for diabetics

Accessories UVL-36, UVL-36 S



Art. 330900 for UVL-36

UV-lamp tubes 1 piece
(9 W, $\lambda = 365$ nm, socket = G23)



Art. 330901 for UVL-36

UV-lamp tubes set of 4 pieces
(9 W, $\lambda = 365$ nm, socket = G23)



Art. 330904 for UVL-36 S

UV-lamp tubes 1 piece
(9 W, $\lambda = 365$ nm, socket = G23)



Art. 330905 for UVL-36 S

UV-lamp tubes set of 4 pieces
(9 W, $\lambda = 365$ nm, socket = G23)

Accessories Nailfan



Art. 230800

filter unit - washable



Art. 230803

arm rest



Art. 230804

mesh grid

Accessories Nailfan mini



Art. 295001

filter unit (How to clean see instruction manual Nailfan mini)

Accessories AC-4000



Art. 334001

filter set



Art. 334002

UVC lamp

Accessories Therapy



Art. 270960

Swiss Depil - replacement pads

Replacement friction pads (6 pieces)



Art. 402091

Arm cuff for PBM-3.5, size L (32-42 cm)

Arm cuff, size L (32-42 cm), grey, for blood pressure monitor



Art. 414003

Arm cuff for BDP-200, size L (33 - 43cm)

Arm cuff, size L (33 - 43 cm), black, for blood pressure monitor



Art. 402093

AC-Adaptor

6 Volt AC Adaptor for blood pressure monitor **PBM-3.2, PBW-5.2** and **BDS-700, PBM-4.2**



Art. 414013

Arm cuff for BDS-700, size M (30 - 42 cm)

Arm cuff, size M (30 - 42 cm), grey, for blood pressure monitor



Art. 414014

Arm cuff for BDS-700, size L (42 - 48 cm)

Arm cuff, size M (42 - 48 cm), grey, for blood pressure monitor



Art. 402061

Arm cuff for PBM-4.2, size (22 - 42 cm)

NEW

Accessories Inhalator



Art. 422001

Protective cap for INH-2.1



Art. 422002

Medicament holder INH-2.1

Spare medicament holder for reordering



Art. 422003

Mouthpiece for INH-2.1

100 % latex-free



Art. 422004

Mask, size S for INH-2.1

suitable for children, with connecting piece

100 % latex-free



Art. 422005

Mask, size L for INH-2.1

suitable for adults, with connecting piece

100 % latex-free



Art. 422006

Adapter for INH-2.1 (picture is given as an example only)



Art. 359174

Neck Electrode

self-adhesive, 1 piece/ pack.; approx. 15 x 9,5 cm

NEW



Art. 359173

Electrode Butterfly

self-adhesive, 2 pcs/ pack.; 15 x 9,5 cm



Art. 359170

Textile electrode

self-adhesive, 4 pcs/ pack; 40 x 40 mm



Art. 359171

Textile electrode

self-adhesive, 4 pcs/ pack; 45 x 45 mm



Art. 359172

Textile electrode

self-adhesive, 4 pcs/ pack; 45 x 35 mm



Art. 359178

Textile electrode

self-adhesive, 4 pcs/ pack; 45 x 80 mm



Art. 359177

Textile electrode

self-adhesive, 4 pcs/ pack; Ø 32 mm



Art. 359070

Carbon permanent electrode

4 pcs./pack + 4 pcs. self-adhesive gel pads;
40 x 40 mm (reorder of individual gel pads, item no. 359280)



Art. 359280

Self-adhesive gel pad

For replacing the used gel pads for the carbon permanent electrodes (item no. 359070)
40 x 40 mm, 10 pcs./pack, self-adhesive

Accessories TENS/EMS



Art. 359910

TENS contact gel

to refresh the contact surfaces of the electrodes. 50 g tube



Art. 359013

Contact cable

contact cable for **EMT-4**, pairs



Art. 359014

Contact cable

contact cable for **EMT-6, T-6 and IT-6**, pairs



Art. 359011

Contact cable

contact cable for **Promed digi 3000, Promed tens 1000 s**, pairs



Art. 359130

Vaginal sonde

Special electrode for vaginal use for the treatment of incontinence and to train the pelvic floor muscles (for women).



Art. 359131

Anal sonde

Special electrode for anal use for the treatment of incontinence and to train the pelvic floor muscles.



Art. 359990

Akku

battery pack/akku for **Promed tens 1000 s, Promed digi tens 3000**, 8,4 V



Art. 359991

Batterie 9 V

battery 9 V for **Promed tens 1000 s, Promed digi tens 3000**



Art. 359995

Akku 4pcs

akku 4 pcs for Promed **EMT-4**
Gr. AAA 1,2 Volt



Art. 359994

Battery 4 pcs

battery 4 pcs for **EMT-4, EMT-6, T-6** Gr. AAA 1,5 Volt

Accessories Blood glucose monitors



Art. M79039

Control solution for MediTouch 2

EAN 40 15588 79039 3



Art. M79038

Test strips for MediTouch 2 (2 x 25 pieces)

EAN 40 15588 79038 6



Art. M79002

Lancing device

EAN 40 15588 79002 7



Art. M79028

Lancets for MediTouch 2 / MediTouch 2 connect

EAN 40 15588 79028 7

Art. M79049

USB cable for MediTouch 2

EAN 40 15588 79049 2

We will gladly send you suitable advertising material:

We constantly renew our range of advertising materials. Contact our employees.

Flyer:

(DIN long - various languages)

(Cardboard) Counter display for DIN long Flyer

Art. 999981



Flyer Electrotherapy (DE)

Art. 993516

Flyer Manicure/Pedicure (DE)

Art. 991218

Flyer Medisana (DE)

Art. M99983

Flyer Shoes (DE)

Art. 995114

Flyer Shoes „Pedibelle“ (DE, EN)

DE: Art. 995401; EN: Art. 995402

Flyer Heating Blankets (DE, EN, FR)

DE: Art. 99977; EN: Art. 994501; FR: Art. 994502

Flyer Personal Care (DE)

Art. 991216

Flyer professional files (DE, EN)

DE: Art. 992020; EN: Art. 992021

Flyer Accessories (DE, EN)

DE: Art. 991901; EN: Art. 991902

Flyer Blood Pressure Monitors (DE)

DE: Art. 994101



Brief instructions (DE):

(DIN A4, cellophane-wrapped)



for TENS/EMS device EMT-4:

Art. 983529



for TENS/EMS device EMT-6:

Art. 983530



Stand-up display for Promed products:

Art. 999901



Promed cloth bag

Art. 999986



Promed ballpoint pen:

Art. 940655



Promed writing pad A6

Art. 999980

A**Abrasive wheel sets:**

Comfort Set	p. 92
Calluses Set	p. 92
Modellage Set Basic	p. 91
Modellage Set HM-Bitset Pro	p. 91
Modellage Set Master	p. 91
Natural Set	p. 92
Personal Care Set	p. 91
Modellage Set Blue-mix	p. 90
Modellage Set Ceramic	S. 90
Modellage Set Rainbow	S. 90

Abrasive wheels:

Blue-Mix Bits	p. 94
Diamond bits	p. 98
Carbide bit, blue series	p. 98
Carbide bit, yellow series	p. 96
Carbide bit, red series	p. 96
Carbide bit, specials	p. 97
Carbide bit, purple series	p. 96
Ceramic bits	p. 93
Polishing bits	p. 101
Rainbow Bits	p. 94
Sapphire bits	p. 99
Steel cutter – stainless	p. 97
Stone and sandpaper bits	p. 99
Titanium bits	p. 95

Air purifier:

AC-4000	p. 06
---------	-------

Aroma diffuser:

AL-150 WS aroma diffuser	p. 45
AL-240 FL aroma diffuser	p. 44
AL-400 CL aroma diffuser	p. 44
AL-1300 WS aroma diffuser	p. 43
AL-305 aroma diffuser	p. 45

Aroma essences/fragrance oils:

Apple-cinnamon (10 ml)	p. 47
Be Happy (10 ml)	p. 46
Cinnamon (10 ml)	p. 46
Cinnamon-orange (10 ml)	p. 47
Eucalyptus (10 ml)	p. 46
Grapefruit (10 ml)	p. 47
Lavender (10 ml)	p. 46
Lemon grass (10 ml)	p. 46
Limette (10 ml)	p. 47
Limone (10 ml)	p. 47
Orange (10 ml)	p. 46
Peach (10 ml)	p. 47
Relax (10 ml)	p. 46
Rose (10 ml)	p. 47
Vanilla (10 ml)	p. 47
Wild cherry (10 ml)	p. 47

B**Blood pressure monitors:**

PBM-3.5 upper arm blood pressure monitor	p. 25
PBM-4.2 upper arm blood pressure monitor	p. 25
HGP-30 wrist blood pressure monitor	p. 26
HGP-40 wrist blood pressure monitor	p. 26

Blood glucose monitor:

Meditouch 2	p. 31
Meditouch 2 connect DUAL	p. 31

C/D**Dust extraction:**

Nailfan	p. 17
Nailfan mini	p. 17

Darco:

AllRound Shoe®	S. 73
AirTraveler Walker	S. 79
Body Armor® Ankle Flex	S. 81
Body Armor® Cast Shoe	S. 76
Body Armor® Night Splint	S. 83
BodyArmor® Heel Reliever	S. 83
BodyArmor® Pro Term	S. 80
Commodus Open®	S. 73
Replacement insole set for WCS®/WCS® Light	S. 83
HeelWedge®	S. 76
MedSurg Pediatric	S. 70
MedSurg Pro® (High)	S. 71
MedSurg Pro® (Low)	S. 71
MedSurg	S. 70
Night Splint®	S. 82
OrthoWedge Light	S. 74
OrthoWedge	S. 74
Posterior Splint Flex	S. 82
Peg Assist Insole®	S. 87
PegContour Insole®	S. 85
PegContour® Insole Uni	S. 86
Puzzle Insole®	S. 84
Relief Contour Insole	S. 84
Relief Dual®	S. 75
Relief Dual® Plus	S. 75
Relief Insert Walker 2.0	S. 79
Relief Insert Walker 2.0 Achill	S. 80
SlimLine® Pediatric	S. 77
SlimLine® overshoe	S. 77
TAS® Toe Alignment Splint	S. 81
Toe Cap	S. 78
TwinShoe®	S. 78
WCS® long-term bandage shoe	S. 72
WCS® Light long-term bandage shoe	S. 72

Disinfection:

Pure-ID (for instruments)	p. 05
UC-50 ultrasound cleaner	p. 16

E/F**Electric files****Professional and Personal Care:**

Emotion	p. 22
Feeling	p. 23
Pedimed	p. 20
Pedisenso Duo	p. 21
Sensitive	p. 22
the file - 620	p. 15
the file - 620 deluxe	p. 20
the file - 625	p. 15
the file - 1030	p. 14
the file - 2020	p. 14
the file - 4030-SX2	p. 13
the file - 4040	p. 10
the file - 5040-SX2	p. 12
Ultra pro S	p. 21

H**Heating products:**

HGP -1.7 heating belt	p. 40
HKP-1.2 H heating pad	p. 41
HKP-1.6 XL, heating pad	p. 41
KFS-100 foot warmer	p. 42
KHP-2.3 cosy heating blanket	p. 38
KHP-2.3 G, heating blanket "Sherpa"	p. 38
NRP-2.0 shoulder and neck heating pad	p. 39
NRP-2.4 heating pad for back and neck	p. 40
NRP-5.4 heating pad for back and neck	p. 39
UBS-2.4 underbed	p. 42

I/J

Inhalation:

- INH-2.1 ultrasound inhalator p. 27

Intimate care:

- Intimate douche "Orchidee" p. 36
- Intimate douche "Rose" p. 36

K/L

Lamps/light:

- LTL-749, table lamp p. 19
- LTM-30, magnifying table lamp p. 19

M/N

Mask:

- MNS-100 mask FFP2 S. 04
- MNS-100 Black, mask FFP2 S. 05

Massage devices:

- ITM Pro hand massager p. 34
- MSA-900 Shiatsu massage seat rest p. 35

O/P

Pulse oximeters:

- PM-200 Fingertip pulse oximeter p. 29
- PM-100 pulse oximeter p. 30
- PM-150 Connect pulse oximeter p. 30

Q/R/S

Shoes:

Pedibelle:

- Pedibelle Alexander p. 56
- Pedibelle Anna p. 55
- Pedibelle Eleine p. 55
- Pedibelle Jonas p. 56

Bandage shoes:

- Flex 2 p. 57
- GentleWalk-Hi p. 61
- GentleWalk-Lo p. 61
- Grainau p. 58
- München 2 LXL p. 62
- München 3 p. 62
- Pantomobil D p. 64
- Sanicabrio DXL p. 59
- Sanicabrio LXL p. 59
- Sanicabrio p. 60
- Sanisoft D p. 60
- Sanisoft p. 61
- TheraLight 1 p. 64
- Theramed D1 p. 63
- Theramed D3 p. 63
- Wallgau p. 58

Soles:

- Replacement soles p. 65
- Huber heel spur p. 67
- Huber joint discharge p. 67
- Huber Hallux Valgus p. 66

Spa:

- FB-100 foot spa p. 35

T

TENS/EMS devices:

- EMT-4 p. 33
- EMT-6 combi unit p. 32
- IT-6 incontinence therapy device p. 32
- Tens 1000 s p. 33

Thermometer:

- IRT-80 infrared thermometer p. 28
- IRT-100 infrared thermometer p. 27
- PFT-3.7 digital clinical thermometer p. 28

U/V

UV lamps:

- UVL-36 S p. 18
- UVL-54 "All-In" p. 18

Z

Accessories:

- AC-4000 accessory p. 105
- Batteries / rechargeable batteries p. 109
- Electrodes p. 108
- Replacement filter bags p. 104
- Replacement discs for Swiss Depil p. 106
- Foot pedals p. 104
- Handpieces p. 103
- Inhaler p. 107
- Contact gel p. 109
- Contact cables p. 109
- Cuffs p. 106
- Nailfan accessory p. 105
- Nailfan Mini accessory p. 105
- Professional accessory kit p. 102
- UVL-36, UVL-36 S p. 105
- Vaginal/anal probe p. 109

- Medisana accessory S. 110

Your contact:

Headquarters/Order processing:

Tel.: +49 (0)8821/9621-0
info@promed.de

Managing Directory, Marketing and Distribution:

Andreas Schuster
Tel.: +49 (0)8821/9621-0

Order processing Export:

Marsha Lo
Tel.: +49 (0)8821/9621-12

Service / Retoure:

Tel.: +49 (0)8821/909 65 50

Business hours:

Monday - Thursday:
8 am - noon and 1 - 5 pm
Friday: 8 am - noon

Address for returns (parcels/returns):

Promed GmbH, Bischofstr. 1
82490 Farchant, Germany

Sales force

Bernd Hegeler

Germany region postal codes:
2, 3, 49
Tel.: 0151/174 815 13

Marlon Osthoff

Germany region postal codes:
4 (without 49); 5 (without 55)
Tel.: 0171/777 71 80

Gabriele Zimmer

Germany region postal codes:
55, 6, 7, 88, 89
Tel.: 0171/932 18 85

Frau Gedak

Germany region postal codes:
80-87; 90-97
Tel.: 0151/113 273 51

Imprint

Changes/GTC

Technical changes further developments, misprints, typographical errors and errors also remain subject to change without notice. You can find our GTC at any time on our website.

Copyright

The information published in the catalogue (text and images, graphics), as well as their structure, is subject to copyright by Promed GmbH, and in parts also the copyright of third parties. Duplication, dissemination, reproduction and transmission, storage and other utilisation is expressly prohibited without our written consent.

© copyright 2023, Promed GmbH, Germany



Promed GmbH

Lindenweg 11
82490 Farchant
Germany

Tel. +49 (0)8821/9621-0
Fax +49 (0)8821/9621-21

E-mail: info@promed.de
Web: www.promed.de

Operation Manager: Manfred Eiberger
Managing Director: Andreas Schuster,
Thomas Hollefeld
Commercial Register B (HRB) No.: 135838
Munich District Court (Amtsgericht)
VAT ID No.: DE 1283/68339

